

Nathanson Liver Retractors ENABLING OPTIMAL ACCESS & FIELD VIEW DURING UPPER GI SURGERY



- Fully supports and retracts the liver for optimal field visualization
- Direct attachment to Mediflex[®] Holding & Positioning Systems
- Hexagonal lock fitting prevents liver slippage and axial rotation
- Reduced trauma compared to hand-held liver retractors
- Extended tip length available for ease of insertion through thicker abdominal walls
- Matte finish to prevent reflection
- Robotic versions to eliminate arm interference
- Widest variety of configurations on the market available in 3mm, 5mm, and 6.5mm diameter and a range of sizes to accommodate all patients

MAXIMUM STABILIZATION • ROBOTIC COMPATIBLE

Nathanson Liver Retractors

Specially designed with hex post fitting for direct attachment to Mediflex[®] Holding & Positioning Systems - for maximum retraction security and stability





RS	Ref. No.	UOM	Dia.(mm)	Size	Tip	Shaft Length (A)
0	69703-R9	Ea	5.0	Small	Std	2" (5.1cm)
C	69704-R9	Ea	5.0	Medium	Std	2" (5.1cm)
IR A	69705-R9	Ea	5.0	Large	Std	2" (5.1cm)
RET	69709-R9	Set	5.0	5.0mm: S, M, L (1 ea)	Std	2" (5.1cm)
<u>0</u>	69736-R9	Ea	6.5	Medium	Std	2" (5.1cm)
OTIC	69736-R12	Ea	6.5	Medium	Std	3.2" (8.1cm)
OB	69736-XL-R9	Ea	6.5	Medium	Extended	2" (5.1cm)
ě	69736-XL-R12	Ea	6.5	Medium	Extended	3.2" (8.1cm)

	Ref. No.	UOM	Dia.(mm)	Size	Tip	Shaft Length (A)
RS	69738-3	Ea	3.0	Extra Small	Std	3.7" (9.4cm)
	69735-3	Ea	3.0	Small	Std	5.5" (14cm)
	69704-3	Ea	3.0	Medium	Std	7.5" (19.1cm)
	69738	Ea	5.0	Extra Small	Std	5.5" (14cm)
2	69703	Ea	5.0	Small	Std	5.5" (14cm)
RETRAC	69704	Ea	5.0	Medium	Std	7.5" (19.1cm)
ETR,	69705	Ea	5.0	Large	Std	10.7" (27.2cm)
	69709	Set	5.0	5.0mm: S, M, L (1 ea)	Std	See above
ARD	69735	Ea	6.5	Small	Std	5.5" (14cm)
DA	69735-XL	Ea	6.5	Small	Extended	5.5" (14cm)
AND,	69736	Ea	6.5	Medium	Std	7.7" (19.6cm)
ST	69737	Ea	6.5	Large	Std	10.7" (27.2cm)
	69710	Set	6.5	6.5mm: S,M,L (1 ea)	Std	See above
	69736-XL	Ea	6.5	Medium	Extended	7.7" (19.6cm)
	69737-XL	Ea	6.5	Large	Extended	10.7" (27.2cm)
	69710-XL	Set	6.5	6.5mm: S, M/XL, L/XL (1ea)	Extended	See above



www.mediflex.com

USA 800.764.0474 +1.631.582.6424 info@mediflex.com

© 2019 Mediflex[®] Surgical Products, Islandia, NY. All rights reserved. Mediflex[®] is a registered trademark of Mediflex[®] Surgical Products, a division of Flexbar Machine Corp.

NLR-0819