

axiostat[®]


stops bleeding, instantly

100%
chitosan
HAEMOSTATIC
dressing

dedicated
to **stop**
bleeding



Axiostat® is a sterile, single use, non-absorbable **hemostatic** dressing intended to be used for temporary control of moderate to severe **bleeding wounds**.




Axiostat® is a 100% chitosan dressing made from pure grade chitosan extracted from shellfish. Chitosan is a linear polysaccharide consisting of glucosamine and N-acetyl glucosamine chains and is highly bio-active. Cationic nature of chitosan makes it mucoadhesive and thereby stops bleeding instantly.

Axiostat is manufactured in **ISO 13485:2016** and **GMP** certified facility.


- Effective on patients on **blood thinners**
- **Active** on all sides
- **Hypoallergenic** no known side-effects
- **Can be cut, folded,** stuffed into deep wounds
- Gamma **Sterilized**
- **Biocompatible**
- Easy to **Apply & Remove**

Axiostat is **India's first** **USFDA** approved Haemostat

HOW IT WORKS



CLINICAL DATA



For detailed Cases/Clinical Studies, Visit:
www.axiobio.com/case-studies or
Scan this Code

AXIOSTAT APPLICATIONS



SURGICAL

L88 - 8x8cm | M85 - 8x5cm
S55 - 5x5cm | V55 - 5x5cm



All ectomy cases



General Surgery

*Suture line anastomosis | Arterial and venous spurting
Cardiothoracic & Vascular Surgeries (CTVS) | GI Surgery
Post resection tumour bed bleeding | Capillary oozing
Donor site bleeding | Orthopedic & Traumatic Surgeries*



VASCULAR

V55 - 5x5cm
V35 - 3.5x3.5cm



Interventional Radiology procedures



Interventional Cardiology procedures (Femoral & Radial)

Any type of Vascular interventional procedures



MILITARY

MIL88 - 8x8cm



Gunshot wounds




Puncture or Stab wounds

Blast injuries | Scalp and Limb cuts | Avulsions



HAEMODIALYSIS

N22 - 2x2cm



CVC procedure



Haemodialysis

*Vascular access procedures | Intravenous catheters
Arterial/Central Line Removal*



EMERGENCY

L88 - 8x8cm | M85 - 8x5cm
S55 - 5x5cm | N14 - 1x4cm
N18 - 1x8cm



Head & Neck Trauma



Arterial bleeding

*Venous Bleeding | Scalp Injuries | Cuts & Lacerations
Limb injuries*



DENTAL

D11 - 1x1cm



Maxillofacial trauma



Oral Surgeries

Dental extractions and other oral procedures

see
bleeding
use **axiostat**[®]

ORDERING INFORMATION

PRODUCT CODE	PRODUCT SIZE	PACKAGING
L88	8cmx8cm	Box of 10
M85	8cmx5cm	Box of 10
S55	5cmx5cm	Box of 10
N18	1cmx8cm	Box of 20, 40
N14	1cmx4cm	Box of 40, 80
MIL88	8cmx8cm	Box of 10
V55	5cmx5cm	Box of 10
V35	3.5cmx3.5cm	Box of 15
N22	2cmx2cm	Box of 20, 80
D11	1cmx1cm	Box of 20, 80



AXIO BIOSOLUTIONS PRIVATE LIMITED

Plot 18, Gujarat Pharma Techno Park, Sari Matoda, Sanand Taluka, Ahmedabad - 382220, Gujarat, India

Customer care: +91 8860786067, Email: info@axiobio.com | www.axiobio.com

twitter.com/axiobio facebook.com/axiobio

www.axiostat.com