


FUJINON FIBEROPTIC BRONCHOSCOPES Specifications



	FB-120S	FB-120T	FB-120P
Viewing direction	Forward		
Observation range	1 ~ 50mm		
Field of view	120 °		100 °
Distal end diameter	4.8mm	5.9mm	2.7mm
Flexible portion diameter	4.9mm	5.9mm	2.8mm
Bending capability : UP	180 °		
DOWN	130 °		
Forceps channel diameter	2.2mm	2.8mm	1.2mm
Working length	600mm		
Total length	900mm		920mm



Light source

PS2-HP

The PS2-HP is a sturdy, reliable lightsource designed for complete portability. A single, 150 watt halogen bulb provides ample illumination for efficient and effective examination. Light intensity can be manually controlled to adjust to all clinical requirements. Additionally, the air pump can be switched on and off in the event that insufflation is not required.

Specifications

Power supply	AC117V/230V
Frequency	50/60Hz
Bulb used	15V, 150W halogen bulb, gold mirror
Dimensions	195 × 350 × 215mm
Weight	8.6kg

FUJIFILM GROUP

FUJINON CORPORATION

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN
 TEL : 81-48-668-2153 FAX : 81-48-668-1570
 E-mail:sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.
 TEL : 1-800-490-0661 FAX : 1-973-633-8818
 E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY
 TEL : 49-2154-924-0 FAX : 49-2154-924-290
 E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,
 SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911
 E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA
 TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD,79 RI JING ROAD,
 WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA
 TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)
 FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com