


max. Ø

 harte Drähte
 hard wires
 alambres duros
 fili duri

max. Ø

 weiche Drähte
 soft wires
 alambres blandos
 fili molli

Tc
 TUNGSTEN

40-876-23 23,5 cm, 9 ¼ "

 1,6 mm 2,0 mm
40-877-16
40-877-16 18,0 cm, 7 "

 1,6 mm 2,4 mm
40-878-18

 längsdurchbohrt
 cannulated
 barrenado oblongo
 tubolare

40-878-18 18,0 cm, 7 "

 1,0 mm 1,5 mm
40-880-14
Tc
 TUNGSTEN

40-880-14 14,5 cm, 5 ¾ "

 1,5 mm 2,0 mm
40-882-14
Tc
 TUNGSTEN

40-882-14 14,5 cm, 5 ¾ "