

# Consona N7

Diagnostic Ultrasound System


# Concrete to Resonate











# Confident

### **Excellent Image Platform**

The ZST<sup>+</sup> platform is an extraordinary innovation, representing an ultrasound evolution. Transforming ultrasound metrics from conventional beam-forming to channel data based processing. It delivers exceptional image quality for infinite imaging solutions with non-stop improvements.






### **Image Gallery**







Supraspinatus Tendon


Kidney CDI with HR Flow

CCA Duplex





Fetal Spine 4D



A2C CDI

Uterine Fibroids




Dedicated to deliver more substantial value in each specific application scenario

## **General Imaging Solution**

Smart tools for smooth diagnosis

#### **Smart HRI**

Providing an easy and fast way of liver steatosis assessment, the higher the HRI value, the deeper the indication of steatosis



Smart HRI




**iTouch**<sup>+</sup> Real-time dynamic image optimization



Smart Calc

Automatic trace, measurement and calculation tool




NTE
ROI stiffness display with 'Shell'
quantitative assessment

### **Women Health Solution**

Full-stack smartness solution based on Smart Scene 3D

#### **Smart Scene 3D**

Full-stack smartness obstetric solution empowered by ZST<sup>+</sup>






Smart Face
Automatic removal of occlusions
over fetal face with one click



Smart OB

Automatic measurement of fetal biometrics




Smart FLC 2D / 3D
Automatic follicle counting and measurement in 2D/3D mode

## **Cardiovascular Solution**

Effective cardiovascular evaluation



AutoEF
Automatic measurement of Ejection Fraction by diastole/systole frames



TTQA

Myocardial synchronization evaluation with quantitative analysis



**RIMT** RF-data based real-time IMT measurement

### **Ergonomic Design**

Boosts the daily work efficiency. Exquisite machine design brings users better operating experience



# Conducive

## Inspires the growth by every progressive step

Shortened learning curve with effective professional improvement

### iScanHelper

A scanning tool to practice skills by following the guidance step by step



### iWorks

A standardized protocol, scanning step by step without plane omission



#### **Smart Vue**

Automatic recognition on standard plane, and restore the standard plane automatically



