

Description of the product

1. The exact name of the product

Treatment table; treatment table with bowls



2. Destiny of the product

Treatment tables may be dedicated to clinics, hospitals and treatment rooms

3. Dimensions and technical datas

Height : 850 mm
Length : 750 mm
Width : 450 mm
Permissible load for shelf : 10 kg

4. The standards met by the product

Declaration of conformity CE

5. Set of equipment fitted

Powder coated frame, 2 shelves made of glass, 4 wheels including 2 with breaks, bowls made of stainless steel

6. Warranty

24 months

7. Producer

UBM – Katarzyna Kotasińska
46-100 Namysłów
ul. Kazimierza Pułaskiego 4B,
NIP 752-111-59-85