



EN-N5 Series Enteral Feeding Pump

EN-N5 Series Enteral Feeding Pump

Basic Parameters	
Case Material	PBT+PC
Screen	4.3 inch LCD color touch screen
Feeding Mode	Continuous, Intermittent, AntiBolus
Accuracy	±5%
Feed Rate Range	1-1200ml/h (min. increment 1ml/h)
VTBI	0-9999 ml
KTO	1-30ml/h adjustable (min. increment 1ml/h)
Bolus Rate	1-1200ml/h (automatic & manual bolus available)
Purge Rate	1-2000ml/h (automatic & manual purge available)
Wash Rate	1-2000ml/h (automatic & manual wash available)
Anti-pumping Speed	1-1200ml/h (min. increment 1ml/h)
Occlusion Level	3 levels selectable (150, 300, 600mmHg)
Interface	DB15 multi-functional interface
Wireless	WiFi (optional)
Heater	30.0-40.0°C adjustable (only suitable for EN-N5H)
Alarm Type	Feeding complete, Pressure high, Battery empty, KTO finished, Dropper empty bottle, Standby time expired, Door Open, Anti-Bolus finish, Temperature high, Temperature low, Nutrient leakage, Heater failure, Nearing completion, Battery near empty, Reminder alarm, No power supply, No heater installed, System error, etc

Special Features	
Titration	Change feed rate without stopping feeding
Reset occlusion level	Reset occlusion alarm level without stopping feeding
Last Therapy	Last therapies can be stored and used for rapid infusion
Anti-Bolus	Automatic drop line pressure to reduce bolus impact after occlusion
Purge	Remove the air bubble
Wash	Flush the nutrition line
Screen Lock	Avoid misoperation
Standby Mode	Standby time adjustable
DPS	Real-time pressure display in graphically and numerically
Drop Detection	Support
Brightness	Adjustable
Alarm Volume	Adjustable
History Record	More than 5000 entries

Power Source	
AC Power	100V-240V, 50/60Hz
External DC Power	10-16V
Battery	Lithium Polymer battery 2600mAh*2, more than 10 hours operating time at 25ml/h

Others	
Classification	Class I, CF
Standard	CE; EN1789:2014 (Ambulance)
Waterproof	IP44
Dimension	139(W)*140(H)*91(D) mm (without Pole clamp)
Weight	Approx. 1.3kg (Pole clamp included)



Shenzhen ENMIND Technology Co., Ltd.

Add: 5th Floor, Block A, Defengsheng
Building, No. 41 Dabao Road, Bao' an
District 23, Shenzhen 518101, P.R. China

Web: www.enmind-tech.com
Email: sales@enmind-tech.com
Tel: 0086 755 2331 6007



Facebook

