



ASU-I

Electric Suction Unit

Features

- Designed for removing fluids obstructions from the airways to facilitate breathing
- Portable electric vacuum pump
- High flow & low noise
- Air filters, over-flow protection
- Oil-free vacuum membrane pump
- Continuous negative pressure regulation
- Ultra-clear vacuum gauge display
- Long-lasting battery
- IP rating: IP43



CE 0051 RTCA/DO-160G EN1789

Parameters

L * W * H	400*240*150mm	Battery capacity	14.8V, 2600mAh
Weight	≤ 5.0kg	Flow rate	24L/min
Operating temperature	+5°C ~ +40°C	Vacuum level	> 80kPa
Storage temperature	-20°C ~ +60°C	Vacuum level accuracy	±5kPa
Operation (storage) humidity	Less than 95%, non-condensing	Collection jar	1L
Internal power supply	DC12V,5A	Noise value	≤ 55dB
External power supply	100-240V~ 50/60Hz	Battery operating time	≥ 40 minutes

Amoul® Ambulanc (Shenzhen) Tech. Co., Ltd.



Building #8, Skyworth Innovation Valley Phase II, TangTou Road No.1, ShiYan, Bao'An District, ShenZhen, China
 +86 755-26072286
www.amoulmed.com
info@amoulmed.com



Amoulmed



ambulshenzhen.Tech



Amoulmed



Amoulmed



Amoul