



# FFsonic UF-750XT

---

Portable Color Doppler  
Ultrasound System  
FFsonic UF-750XT

---

## Features

- Fully software controlled Image Processing (F-XT technology)
- Color Doppler & Power Doppler display available
- Autofocus function for transmission
- High-resolution SVGA display
- Cine memory with 128 frames
- Wide band Probe for multi frequency function
- 6-step STC adjustment
- DICOM-compatible Ethernet port
- 3.5-inch 640MB drive for image filing

---

### Dimensions:

Approx. 26(W) x 18.2(D) x 6.05(H) cm  
Approx. 15(W) x 8 1/2 (D) x 14 1/2(H)  
inch

### Weight:

Approx. 2 Kg  
Approx. 26 1/2 lbs

---



**FFsonic**  
**UF-750XT**

