

Full Digital Ultrasound System UF-850XTD



Features

- Full Digital Beam Former
- Advanced F-XTD Imaging Technology
- Multi Beam Processing
- High Density and Wide Band Probes
- Quick Boot up Time
- Wide Flat Panel Display
- Active 3 probe ports
- Thin and flexible probe cable within pin-less connector
- CWD mode Display (option)
- Triplex Display mode
- High frequency Imaging up to 14MHz
- User-friendly Interface with backlit intuitive keys
- Multiple Tissue Harmonic Imaging
- Digital Image workstation provided
- Digital Motion image Recording capability (option)

Applications

- Abdominal
- OB/GYN
- Cardiac
- Vascular
- Musculoskeletal
- Small Parts
- Pediatrics
- Neonatal
- Urology
- Veterinary
- Intraoperative
- Biopsy

Dimensions:

Approx. 48(W) x 79 (D) x 136 - 145 (H) cm

Approx. 19 (W) x 31 (D) x 53 1/2 - 57 (H) inch

Weight:

Approx. 100 Kg

Approx. 220 1/2 lbs