

## Equipment profile

Manufacture	ESAOTE
Name of the equipment	Portable Ultrasound
Type	MyLab Touch
Origin	Switzerland
Equipment ID	506
Manufacturing date	2010
State of quality	Demo device, Excellent, as is
Numbers of equipments	1

### Description

#### MyLab One

The MyLab One is designed to bring ultrasound at the point-of-care. It represents your arm-held solution: steadily fix the ultrasound device to your arm, providing comfortable and stable setting, and it will be able to face any working condition, even the most critical ones.

Emergency and critical care can require very high standards of technology: time and reliability of diagnosis are crucial factors. MyLab One integrates all these characteristics ensuring a long working time with battery. It offers the highest standard of ergonomics, both for the system itself, and for the wide range of accessories and probes.

#### Features:

- Innovative 12" full touch-screen high-resolution TFT LCD color monitor
- Complete multi-frequency probes range
- Superior image quality with unmatched high-frequency range
- PW Doppler with automatic Doppler tracing
- Extreme portability and flexibility of use
- Portrait-landscape auto-orientation, controls on the handle, belt and wrist locker
- TEI™ - Tissue Enhancement Imaging
- Tp-View for enlarged field of view, specially in high-frequency and vascular exams
- MView - combined standard and steered ultrasound imaging to detect all anatomical structures
- XView - eXtreme View with real-time adaptive algorithm
- RF-QIMT & RF-QAS - Quality Intima Media Thickness for early detection of cardio-vascular diseases & Quality Arterial Stiffness for early detection of cardio-vascular diseases
- NNE, Nerve & Needle Enhancement - dedicated Anaesthesiology technology & presets
- Advanced Cardiac Package
- CFM/Power Doppler capability
- Fast start-up time, long life integrated battery, battery recharge, additional battery package



#### Nasz adres w Polsce:

MedTechTrade AG  
ul. Złota 61  
00-819 Warszawa  
Phone +41 43 399 01 11  
Fax: +41 44 905 19 29  
contact-poland@medtechtrade.com



#### Office Schweiz:

MedTechTrade AG  
Bankstrasse 4  
CH-8610 Uster  
Switzerland  
Phone +41 43 399 01 11  
Fax: +41 44 905 19 29  
contact@medtechtrade.com



#### Oficinas Perú:

MedTechTrade Perú S.A.C.  
Calle Juan Fanning 219  
Miraflores (Lima)  
Perú  
Tel: +51-1-447-2115  
contact-peru@medtechtrade.com



#### Oficinas Chile:

MedTechTrade Chile  
Monte Carmelo 37  
Providencia  
Santiago, Chile  
Fono: +56-2 710 5560  
contact-chile@medtechtrade.com

- Wireless networking for easy data sharing, storage and printing
- Remote archive, review, report, back-up
- Up-to-date connectivity and data management solutions (wireless, LAN, integrated database)

Incl. Trolley

Device for General Imaging use.

Probes: Convex

Price EXW Uster/Switzerland

For any question do not hesitate to contact us. We will assist you in German, English, Spanish, Italian, French, Polish

\*Picture for reference only



**Nasz adres w Polsce:**

MedTechTrade AG  
ul. Złota 61  
00-819 Warszawa  
Phone +41 43 399 01 11  
Fax: +41 44 905 19 29  
contact-poland@medtechtrade.com



**Office Schweiz:**

MedTechTrade AG  
Bankstrasse 4  
CH-8610 Uster  
Switzerland  
Phone +41 43 399 01 11  
Fax: +41 44 905 19 29  
contact@medtechtrade.com



**Oficinas Perú:**

MedTechTrade Perú S.A.C.  
Calle Juan Fanning 219  
Miraflores (Lima)  
Perú  
Tel: +51-1-447-2115  
contact-peru@medtechtrade.com



**Oficinas Chile:**

MedTechTrade Chile  
Monte Carmelo 37  
Providencia  
Santiago, Chile  
Fono: +56-2 710 5560  
contact-chile@medtechtrade.com