

- Small & Light Weight Design
- Digital LEDs Display
- Double-Direction Display Mode
- Displays: SpO₂, Heart Pulse Rate and Heart Frequency Bar
- Low Power Consumption (30 Hours Non-Stop Operation)
- Low-Battery Indicator
- Automatic Shut-Off
- Suitable for Several Finger Sizes (from Pediatric to Adult)



Available from January 2011

STANDARD ACCESSORIES

- 2 x AAA 1.5V Alkaline Batteries (Replaceable)
- 1 x Hang Lace

TECHNICAL FEATURES

Patient Range	Adult and Pediatrics Patients
SpO ₂	Measurement range 70% - 99% Resolution 1% Measurement accuracy 80% - 99%: ±2%; 70% - 79%: ±3% ≤69%: no definition
Pulse Rate	Measurement range 30 - 235 bpm Resolution 1bpm Measurement accuracy ±2 bpm or ±2%
Size	cm 5,8 x 3,2 x 3,4
Weight	50g (including two AAA batteries)
Power Source	2 x Alkaline Batteries - AAA size, 1.5V
Warranty	2 Years
Units per Shipping Carton	40
Ref. Number	OX 100000

CE 0123

