

NEW SPEEDYMED

NEW SPEEDYMED is a piston compressor nebulizer for aerosoltherapy. Designed for frequent treatments in home-care, clinics or hospitals with large compartment for convenient storage of accessories and medications.



Lite ref.:

RE 300100

STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
- Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)
- Air Tube
- Air Filter

For spare consumables see page 27

CE 0123

TECHNICAL FEATURES

Compressor	Oiless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
Power consumption	184 VA
Max Pressure	2.5 bar 36 psi 250 kPa
Operating Pressure	1.30 bar 18,9 psi 130 kPa
Max Air Flow	16 l/min
Operating Air Flow	5,2 l/min
Neb-rate with 4 ml solution	0,50 ml/min
Mass Median Aerodynamic Diameter	3,25 µm (MMAD) according to EN-13544-1
Respirable Fraction (< 5µm)	> 70%
Duty Cycle	Non-stop operation
Noise Level (EN 13544-1 Standard)	Approx. 57 dBA
Weight	2,25 kg
Size (Main Unit)	cm 23 x 21 x 15
Years of Warranty	2
Place of Manufacturing	Italy
Shipping carton	6