



ACP-31

Auto Chart Projector





The newly upgraded and multi-functional Auto Chart Projector ACP-31 by our advanced technology.


DONGYANG's New Auto Chart Projector ACP-31 is our successful achievement for refraction examination with maximum care for high resolution image.

New Auto Chart Projector **ACP-31** Provides you with the function of which remote controller is available to select one of the signal channels in order to prevent other signals from receiving.


SPECIFICATIONS

Charts	31 pieces	
Masks	1 Open / 5 Horizontal Lines / 5 Vertical Lines / 21 Single Letters	
Chart Selection Speed	Average 0.2sec	
Filter	Red / Green	
Projection	Distance	1.5~10.0m (5m is standard)
	Magnification	30× (at 5m)
	Lamp	12V 35W Halogen
Tilting Angle	Ball Joint (±12degree)	
Power Supply	AC 110V, 200-240V 50/60Hz	
Power Consumption	50W	
Power Saving Function	5 Minutes After	
Interface	9 Pin Connector for Simultaneous Control of Chart	
	Mask and R/G filter by a motorized Refractor's Control unit	
Accessories	Dust Cover, Polarized Metal Screen, Fixation Holder, Polarized Glasses,	
	Spare Fuses, Batteries, Spare Lamp	
Dimension	260×230×360mm	
Weight	Approx. 4.0kg	
Remote Controller	Use of 5 selecting Channels	

Optotypes



A type



B type



Remote Controller



DONGYANG Optics Co., Ltd.

#41-7 Wonmi-dong Wonmi-gu Buchon-city Kyunggi-do Korea Tel. +82-32-613-4317
Fax. +82-32-613-4319 Email: trade@dyopt.com Website <http://www.dyopt.com>