



SL-R5Z Slit Lamp

The SL-R5Z is one of Ray Vision' s NEW 3rd Generation Slit Lamps. Combined traditional design with superior optics, SL-R5Z will be a welcome addition to your practice.

- Classical Zeiss type slit lamp
- Outstanding Galilean optics with moistureproof, mildewproof and anti-reflection lens
- 0 mm to 14mm slit width
- Drum change 6X, 10X, 16X , 25X and 40X magnifications
- Optional external illumination system
- Single-handed elevation system

SL-R3Z Slit Lamp

The SL-R3Z is one of Ray Vision' s NEW 3rd Generation Slit Lamps. It can be used as a conventional examination instrument with excellent performance and high quality.

- Classical Zeiss type slit lamp
- Outstanding Galilean optics for exceptional clarity and detail
- 0 mm to 14mm slit width
- Three magnification ratios: 10X, 16X, 25X
- Versatility and ease of operation

