






eZono™ 3000

Ultrasound for Anaesthesia & Emergency



- Portable & robust
- B and M-Mode
Color & Power Doppler
- Image & Video
storage to USB
- Easy & sterile use
via touch screen



Patented Cue Cards

- Standardized workflow and image interpretation
- Accelerate the learning curve on nerve blocks & vascular access
- Enhance patient education

Leading Value Made in Germany




Image Quality

Exclusive to the eZono 3000™, the user can choose the image style they prefer



eZono™ 3000 Probes

Range of linear, convex and micro-convex probes available with unique features for:

Nerve Blocks e.g.

- Plexus Brachialis
- Sciatic Nerve
- Epidural & Spinal
- Rescue Blocks

Emergency Procedures e.g.

- FAST
- FEEL
- PREP
- Abdominal Aortic Aneurysm
- Pleural Examination



eZono™ 3000 Specifications

System:

- Weight: 4,4kg
- Dimensions: 32,4cm x 28,2cm x 9,3cm
- 4GB internal memory
- 26,4cm diagonal LCD touch screen (24bit, 60 million colors)
- Digital broadband architecture
- 2x USB 2.0 ports, 1 Express Card port, built-in microphone.

Software:

- Proprietary operating system
- less than 10 secs startup time for instant availability
- Icon-based graphical user interface via touch screen
- Upgrade via USB for continuous improvement throughout its life time
- Store up to 30,00 images or 30 minutes of continuous video

eZono AG
Spitzweidenweg 30
D - 07743 Jena

Telefon: +49 (0)3641 876 1730
Telefax: +49 (0)3641 876 1759
sales@eZono.com
www.eZono.com

With compliments

eZono

empowering Sonography®

Document Version 1.0