

MeriAqua Shower Trolley

Safety and ease of use when showering and transferring patients



- **light to maneuver**
- **safe functions**
- **patient comfort**
- **hygienic, high quality materials**
- **long service life**

The MeriAqua shower trolley is the perfect piece of equipment for patient showers and transfers. The shower mattress is soft against the skin and the folding sides make patient transfer easy. The grooves on the mattress surface carry out the wash-water to a drain opening and through the tubing to the floor drain. The tilt adjustment allows water to drain quicker. The MeriAqua shower trolley is also equipped with height adjustment, which allows the nurse to adjust the basin height to an ergonomically suitable level. The central braking system and easy-glide castors make the trolley easy to maneuver. All the handles and levers used to operate the trolley are located at convenient and safe positions within an easy reach.

MeriAqua Shower Trolley



Showermattresses, drainage hose and cushion are available as an accessory.



The side of the mattress folds down, which makes patient transfers easier. The side rails lock easily and safely when raised up. The mattress material holds up well to frequent cleaning.



The tilt adjustment allows for more efficient water drainage.



The clear design of the shower trolley makes it easy to clean.



The tilt adjustment lever is within easy reach. The drainage hose detaches easily for cleaning.



The easy-glide specialty casters and the convenient central brake system provide extra safety to the shower trolley.

Technical Specifications

- Frame of stainless steel
- Mattress base of antistatic ABS-plastic, size 660 x 1955 mm
- Swing-down side rails of stainless steel, lock automatically in place when raised
- Castors Ø 150 mm, equipped with central braking system
- Max. load 150 kg

Dimensions

- Max. width 780 mm
- Max. length 2050 mm

Adjustments

- Height adjustment hydraulically 600 - 900 mm, pedals on both sides of the trolley
- Mattress base longitudinal tilt adjustment (4°) stepless, gas spring aided

Accessories

- shower mattress, special plastics, padded
- cushion
- drainage hose



Merivaara Corp.
Puustellintie 2, FI-15150 LAHTI, FINLAND
Tel. +358 3 3394 611 • Fax +358 3 3394 6144
merivaara@merivaara.fi

www.merivaara.com