

Ada

The most easy-to-use reading and dinner table adjuster mechanism on the market:

- can be adjusted with one hand
- stepless height adjustment
- handedness can be easily and quickly changed without any tools



The twin castors make the cabinet easy to move. Castor locks anchor the cabinet firmly in place, a drawer lock is standard equipment. The cabinet is also available without the reading/dinner table (Figure a). With a wide range of accessories, the cabinet can be modified to meet specific needs.

The Ada Bedside Cabinet frame is made of epoxy-coated steel, top and side panels are light grey, birch or beech melamine. There are three colour options for the drawer fronts. The plastic drawers are lightweight and easy to clean.

Technical specifications and accessories

Dimensions

- width: 475 mm (dinner table width: 545 mm)
- height: 900 mm
- depth: 600–1130 mm (with dinner table in front)
- castor Ø: 75 mm (individually-locking)

Side panels

- grey melamine.....10000 9186
- birch melamine10000 9187
- beech melamine.....10000 9188

Drawers

- white fronts10000 9210
- light blue fronts10000 9212
- blue fronts.....10000 9218

Dinner tables

- reading and dinner table.....10000 9190

Accessories

- table top cover ABS plastic (b)10000 9262
- magazine basket (f).....10000 9266
- slipper rack (b,d)10000 9268
- patient intercom unit holder10000 0629
- bottle holder for lower drawer (e).....10000 9270
- tray for small items, narrow (e).....10000 9280
- tray for small items, wide (f).....10000 9285
- bottle or waste bag rack (g).....10000 9288

warranty 3 years

