

Nova



2670



2660

Nova Bedside Cabinet 2670 and 2660 have easy-to-clean surfaces. Twin-castors ensure easy and quiet movability.



2660

Nova Bedside Cabinet including reading/dinner table. The balanced construction makes the bedside cabinet safe when using the reading/dinner table. Adjustments can be conveniently made regardless of the bed's position.



2781

Separate holder for patient intercom unit can be attached with finger screw. Colours according to the drawers.



2680

Nova Reading/Dinner Table 2680. The Nova collection includes a separate stand-mounted dinner/reading table.

Nova – a stylish bedside cabinets are functional and efficient to hospital and nursing home care facilities. Modern materials and manufacturing technology make the bedside cabinets easy to handle and keep clean. The sturdy, adjustable reading/dinner table can be accessed from either side of the cabinet. The easy gliding castors ensure that the bedside cabinet moves easily and quietly. Rounded edges and reliable functions ensure safe use.

Technical specifications and accessories

- Frame of anodized aluminium profile and composite material
- Drawers made of melamine or wooden laminate (max. load 50 kg)
- Twin castors Ø 50 mm/75 mm, separate castor locks
- The height of reading/dinner table steplessly adjustable with gas springs
- Reading/dinner table made of composite material (max. load 10 kg)

Dimensions

- Width: 2660: 660 mm (reading/dinner table lowered)
2670: 480mm
2680: 790 mm/620 mm (reading/dinner table lowered)
- Depth 2660: 880 mm (75mm Ø castors)
2670: 780 mm (50mm Ø castors)
- Height 2660: 880 mm (75mm Ø castors)
2670: 780 mm (50mm Ø castors)
- Reading/dinner tray (model 2660 and 2680)
520 × 360 mm (max. load 10 kg)
- Adjustment range of the reading/dinner table (model 2660 and 2680) 795–1075 mm with Ø 75 mm castors

Accessories and options

- table top cover, ABS-plastic 120002762
- holder for patient intercom unit:
 - stained rose-wood laminate 12000 2794
 - stained beech laminate 12000 2792
 - birch laminate 12000 2790
 - light grey melamine 12000 2781
- cabinet colour options:
 - stained rose-wood laminate
 - stained beech laminate
 - birch laminate
 - light grey melamine

warranty 3 years



A stylish bedside cabinet series