



Examination table 417 and step 446

Examination table 417 is a basic model for general and gynaecological examinations and minor procedures. Table includes stepless adjustments for head- and back section. Adjustable heel supports, paper roll holder and a stainless steel bowl are also included in the price.

- economical model, easy to use
- table includes stepless adjustments for head- and back section
- upholstery easy-to-clean artificial leather, several colour options available (see available colours on the left page)
- includes adjustable heel supports, stainless steel bowl and paper roll holder (paper roll max width 525 mm)

Technical Specifications

Mattress base	3-sectional, epoxy coated steel frame
Length	1895 mm
Width	645 mm
Height	850 mm
Weight	44 kg
SWL	135 kg

Adjustments

Head section	- 90° ... + 48°
Back section	+ 62°

100000417 examination table 417

Accessories

100000428	knee supports 428, pair
100000427	clamp 427 for fixing knee supports 428
100000418	accessory rails 418, pair
100000446	mounting step

Warranty 3 years



Gynaecological position



The knee supports with fixing attachments



provided with specified accessories conforms to the European Council Directive 93/42/EEC, Medical Device Directive, Annex VII and the corresponding Finnish National Law and regulation.

The product is controlled in accordance with ISO 9001 and ISO 13485 Quality Management Systems and ISO 14001 Environmental Management System.