



Merivaara's examination table 408 is an economical, basic model for general examination and small procedures with stepless head section adjustment.

- practical and easy to use
- for general examinations and small operations
- upholstery easy-to-clean artificial leather, several colour options available
- stepless adjustment of the head section angle included
- equipped with a paper roll holder (paper roll max width 525 mm)

Technical specifications

Mattress base	2-sectional, epoxy coated steel frame
Weight	ca 30 kg
SWL	135 kg
Length	1900 mm
Width	655 mm
Height	745 mm







Adjustments

Head section - 20° ... + 74°

100000408	408 examination table, 2-sectional, including paper roll holder+head section adjustment
100000446	mounting step

Warranty 3 years

Examination table colour options:

	Beige 0222
	Green 8549
	Black 0000
	Blue 2819
	Bluish Grey 0887
	Grey 0705



provided with specified accessories conforms to the European Council Directive 93/42/EEC, Medical Device Directive, Annex VII and the corresponding Finnish National Law and regulation.

The product is controlled in accordance with ISO 9001 and ISO 13485 Quality Management Systems and ISO 14001 Environmental Management System.