




# ECG-6010

Digital 6-channel ECG

# ECG-6010

## Digital 6-channel ECG



Optional SD Card



Lan Port



Foldable Handhold



ECG Information System



### Features:

- Digital isolation technology and signal processing ensured the signal reliability;
- 5.7 inches high-brightness LCD;
- Detailed analysis report;
- Economy mode available help to reduce printing cost;
- Build-in rechargeable Li-ion battery;
- Optional SD card, unlimited data storage;
- LAN port for real-time data transmission;
- Multi-language options ;
- Easy upgrading keeps you updated;
- Delightful operating experience with durable silicon button, unique shuttle design.

### Specifications:

Lead	Standard 12 leads simultaneous acquisition	Power Supply	AC:100~240V, 50/60Hz,42VA
Display	5.7 inches LCD		DC:14.8V rechargeable battery
Language Option	Chinese, English, Russian, French, Spanish, Italian, Polish, Czech, Portuguese	Operation Mode	Automatic, manual and rhythm
Input Impedance	>50MΩ	Recording Paper	112mm×27m(roll)
Frequency Response	0.05~150Hz		or 112mm×140mm×20m(fan- fold)
CMRR	>100dB	Safety Level	IEC60601-1:2005 Class I Type CF
Sensitivity	2.5, 5, 10, 20mm/mV	Dimension of packing	400 mm×355 mm×280mm(L× W×H)
Paper Speed	5, 6.25, 10, 12.5, 25, 50mm/s	Gross Weight	6.5 Kg

**Shenzhen Biocare Electronics Co., Ltd.**



Add: Floor 6, Main Building, High-tech Innovation Park Baoan, Shenzhen, 518102, P.R.China

<http://www.biocare.com.cn>

Tel: +86 755 27960641

Fax: +86 755 27960643

E-mail: [sales@biocare.com.cn](mailto:sales@biocare.com.cn)