

FM-200

Fetal Doppler



FM-200 Fetal Doppler

Water-proof Probe



Back light LCD



Operation buttons



Carry bag



Features:

- Handheld & elegant design;
- Back light LCD;
- Display of FHR, work mode, battery status, signal indication and probe frequency;
- Three work modes: real time, average and manual;
- Low power consumption design;
- FHR alarm limits is adjustable and saved;
- Unique heart rate algorithm enable FHR more accurate and reliable;
- Buckle design for easy take;
- Automatically power off while no signal within 2 minutes;
- Inter-changeable water-proof probe of 2MHz or 3MHz.

Options: 3MHz probe, Carry bag, Earphone, rechargeable battery.

Specifications:

Ultrasonic Work Frequency: 2MHz, 3MHz(Optional);
Heart Rate Range: 50bpm~240bpm;
Alarm Limit: Adjustable lower limit range: 50bpm~120 bpm;
Adjustable upper limit range:160bpm~240 bpm;

Power Requirement: 2 x 1.5V (AA) Alkaline Battery
or Rechargeable Battery;

Dimension: 130mm×89mm×34mm;
Net Weight: 385g(including battery).

Shenzhen Biocare Electronics Co., Ltd.

Add: Floor 6, Main Building, High-tech Innovation Park Baoan, Shenzhen, 518102, P.R.China
Tel: +86 755 27960641

Fax: +86 755 27960643

E-mail: sales@biocare.com.cn



<http://www.biocare.com.cn>