



ECG-3030
ECG-3010

Digital 3-channel ECG

ECG-3030

Digital 3-channel ECG

- Features:
- 112mm wide recording, effectively prevent waveforms overlapped;
 - 5.7 inches high-brightness LCD, contrast is adjustable;
 - Detailed analysis report and representative beats;
 - Economy mode available help to reduce printing cost;
 - Optional SD card, unlimited data storage;
 - LAN port for real-time data transmission;
 - Easy upgrading keeps you updated;
 - Delightful operating experience with durable silicon button, unique shuttle design.



Optional SD Card



Lan Port



Foldable Handhold



ECG Information System



ECG-3010

Digital 3-channel ECG

- Features:
- Digital isolation technology and signal processing ensured the signal reliability;
 - High resolution LCD, brightness and contrast are adjustable;
 - Economy mode available help to reduce printing cost;
 - Build-in rechargeable Li-ion battery;
 - Optional SD card, unlimited data storage;
 - USB port for real-time data transmission;
 - Easy upgrading keeps you updated;
 - Elegant design small size and premium silicon button make big difference.



🕒 Detailed Analysis Report



📁 Optional SD Card



👉 Inner Handhold



🖥️ ECG Information System



Specifications:

	ECG-3010	ECG-3030
Display	3.8 inches	5.7 inches
Language Option	Chinese, English, Russian, Spanish	Chinese, English, Russian, French, Spanish, Italian, Polish, Czech, Portuguese
Acquisition Mode	Simultaneous 12-lead acquisition	
Frequency Response	0.05~150Hz	
Patient Leak Current	<10 μ A	
Time Constant	\geq 3.2s	
Noise Level	<15 μ Vp-p	
Input Impedance	>50M Ω	
Sensitivity	2.5, 5, 10, 20mm/mV	
Recording Paper Specification	63mm \times 20m roll type	112mm \times 27m roll type or 112mm \times 140mm \times 20m fan- fold paper
Paper Speed	6.25, 12.5, 25, 50 mm/s	5, 6.25, 10, 12.5, 25, 50 mm/s
Recording Format	Manual / Auto 1, 1+1, 3, rhythm	Manual / Auto. 3, 3+1
Power Supply	AC: 100-240V, 50/60Hz, 60VA DC: sealed maintenance-free rechargeable battery continuous working time \geq 3 hours	
Dimension of packing	360mm \times 310mm \times 225mm	400mm \times 355mm \times 280mm
Gross Weigh	5.0Kg	6.5Kg



ECG-3030



ECG-3010