



Breaking with tradition

ComfortLite 2 gives patients a sensible alternative.

For patients who struggle with compliance using traditional nasal masks, we have a solution to help them be comfortable and compliant.

ComfortLite 2 was inspired by the men and women who have difficulty getting used to conventional nasal masks – or prefer something nontraditional. With this unique design, we offer the choice of three different cushion types – a pillows cushion, the Simple Cushion, and Direct Seal cushion – for a more personal, more comfortable fit.

And, because we use the “halo” headgear and forehead mechanism for stability and adjustment, the mask barely touches the face. This allows the cushion to provide a reliable seal with no pressure points. FitPacks and DuoPacks include a variety of styles and sizes in various combinations, making the mask a good choice for you and your nontraditional patients. ComfortLite 2 is another great example of how Philips Respironics turns ideas and insights into practical solutions that help you fit 100% of your patients, 100% of the time.

PHILIPS

RESPIRONICS

The 360° crown swivels
move easily with the patient.

Multiple position angle adjuster
and rotating vertical adjuster
provide optimal cushion
placement for seal and comfort.

Flexible adjustment band on
pillows cushion and Direct
Seal optimizes cushion angle
for the best seal.



Tubing attachment point sits toward the
forehead, minimizing contact with the
pillow – and maximizing stability.

Soft headgear fabric eliminates pressure
points and improves comfort for all sleepers.

Three nasal interface options



The pillows cushion inserts into
the nostrils for a quick seal and
easy patient set-up.



Simple Cushion is a tip-of-the-nose
nasal cushion that is smaller to
reduce bulk and pressure points.



Direct Seal cushion offers
a comfortable seal without
insertion into the nostrils.

Item	With headgear	HCPCS	Mask only	HCPCS
ComfortLite 2 mask				
Mask subassembly – no cushions	1025225	A7035*	1025194	*
FitPack/DuoPack sizes: Pillows: S, M; Simple: S, M	1030494	A7034, A7035	1030496	A7034
FitPack/DuoPack sizes: Pillows: M, L; Simple: M, L	1030495	A7034, A7035	1030497	A7034
FitPack/DuoPack sizes: Direct Seal: 4, 5; Simple: S, M	1025112	A7034, A7035	1025121	A7034
FitPack/DuoPack sizes: Direct Seal: 5, 6; Simple: M, L	1025113	A7034, A7035	1025122	A7034
FitPack with pillows cushion only; sizes: S, M, L	1030493	A7034, A7035	1030498	A7034
FitPack with Simple Cushion only; sizes: S, M, L	1025111	A7034, A7035	1025120	A7034
FitPack with Direct Seal cushion only; sizes: 4, 5, 6	1025127	A7034, A7035	1025123	A7034
Parts and accessories	Part number	HCPCS	Part number	HCPCS
Halo headgear (including strap)	1025117	A7035	Direct Seal cushion – size 2	1021806 A7033
Halo headgear (strap only)	1025114	N/A	Direct Seal cushion – size 3**	1021807 A7033
Pillows cushion – petite	1030502	A7033	Direct Seal cushion – size 4**	1021815 A7033
Pillows cushion – small	1030503	A7033	Direct Seal cushion – size 5**	1021816 A7033
Pillows cushion – medium	1030504	A7033	Direct Seal cushion – size 6**	1021817 A7033
Pillows cushion – large	1030505	A7033	Simple Cushion elbow (5/pack)	1014206 A7045
Simple Cushion – petite**	1021835	A7032	Pillows cushion sizing gauge (5/pack)	1032007 N/A
Simple Cushion – small**	1021830	A7032	Pillows cushion sizing gauge (20/pack)	1031154 N/A
Simple Cushion – medium**	1021823	A7032	Simple Cushion sizing gauge	1018624 N/A
Simple Cushion – large**	1021821	A7032	Direct Seal cushion sizing gauge (20/pack)	1018665 N/A
Direct Seal cushion – size 1	1021802	A7033	Small accessory pouch	1027999 N/A

*Requires addition of mask cushion for HCPCS A7034.

**Also available in 5/packs.

Respironics, ComfortLite, Direct Seal, DuoPack, FitPack, and Simple Cushion are trademarks of Respironics, Inc. and its affiliates, all rights reserved.



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoehc JJ 09/14/09 MCI 4102595 PN 1028746

Philips Healthcare is part of Royal Philips Electronics

Respironics Asia Pacific
+65 6298 1088
Respironics Australia
+61 (2) 9666 4444
Respironics Europe, Middle East, Africa
+33 1 47 52 30 00
Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics