



## **SERIES-2 BRONCHOSCOPE EFER-DUMON: BRONCHIAL ENDOSCOPES**

**Endoscopes compatible with modular BRONCHOSCOPE EFER-DUMON instrumentation.**

### **Specifications :**

- **3, 4 and 5.5 mm diameters**
- **Axial viewing**
- **32 mm diameter standard eyepiece**
- **Stainless-steel**
- **Removable connectors for compatibility with other brands of light cables**

**Vocation : Medical Endoscopy**

## BRONCHIAL ENDOSCOPES

| REFERENCE      | DIAMETER | WORKING LENGTH | DIRECTION OF VIEW | FIELD OF VIEW |
|----------------|----------|----------------|-------------------|---------------|
| <b>BX-5500</b> | 5.50 mm  | 49 cm          | 0°                | 60°           |
| <b>BX-4000</b> | 4.00 mm  | 49 cm          | 0°                | 60°           |
| <b>BX-3000</b> | 3.00 mm  | 27 cm          | 0°                | 45°           |



### ACCESSORIES

#### Supplied :

C1 light guide connector  
(EFER/WOLF standard)

C2 light guide connector  
(STORZ standard)

#### Optional :

C4 light guide connector  
(Olympus / ACMI standard)

### OTHER PRODUCTS OF THE MEDICAL RANGE



**STANDARD RANGE** - Autoclavable endoscopes, Xenon light sources and digital technology video cameras are the backbone of a smart medical endoscopy line.

### OTHER INSTRUMENTS OF THE SERIES-2 BRONCHOSCOPE EFER-DUMON RANGE

**SERIES-2 BRONCHOSCOPE EFER-DUMON : Adult instrumentation**

**SERIES-2 BRONCHOSCOPE EFER-DUMON : Pediatric and stent placement instrumentation**

**SERIES-2 BRONCHOSCOPE EFER-DUMON : Adult stent placement instrumentation**

EFER reserves the right to modify the product specifications without notice

**EFER**  
ENDOSCOPY