

# COMBO

SLT · Nd:YAG

Two complete laser-systems  
combined in one slit lamp





Nd:YAG Laser



SLT Laser



Two units on one Slit lamp. The SLT and the YAG laser are incorporated into the Haag Streit BQ 900® slit lamp without compromising the quality of the optics.

There is always a place for this Combo - even in the smallest practice.

### SPECIFICATIONS ND:YAG

Laser, wavelength	Nd:YAG
Repetition rate	1 Hz
Pulse length	4 ns
Aiming beam	Laser diode 635 nm, red
Micro offset	50/125/250 µm posterior
Burst mode	2 and 3
Output energy	27 mJ max. at 3 pulses
Energy adjustable in steps	0,5/0,7/1,0/1,3/1,6/2,0/3,0/5,0/7,0/10 mJ in 10 steps
Mains supply	85 to 230 V

### SPECIFICATIONS SLT

Laser, wavelength	Q-Switch Nd:YAG, 532 nm
Focus spot	400 µm
Repetition rate	1 Hz
Pulse length	3 ns
Output energy	0,2 to 2,0 mJ
Energy adjustable in steps	0,2/0,4/0,6/0,8/0,9/1,0/1,1/1,2 /1,5/2,0 mJ in 10 steps

Technical data is subject to change without notice. We always want to maintain the highest standard.

