

Dcare® **medicool 14 P F1**

TECHNICAL DATA



CARE AND READING LIGHTS

FITTED WITH	1 x compact fluorescent lamp TC-TSE 14W/827 E27 - D68880000
LUMINOUS POWER	450 Lux / 0,5 m
COLOR RENDERING INDEX (CRI)	> 80
COLOUR TEMPERATUR	approx. 2700 K
WORK EQUIPMENT	integrated in spare lamp
CONNECTED LOAD	230-240V; 50Hz
MAINS LEAD	3 m Euro
LUMINAIRE BODY	polycarbonate (PC)
WEIGHT (NET)	approx. 1.4 kg
POWER CONSUMPTION	approx. 14 W
BALANCE OF ARTICULATION JOINTS	spring
CLASS OF PROTECTION	II
NORMS	EN 60598-1, EN 60598-2-25
SPECIAL FEATURES	separate light strip for orientation
FASTENING	rail, table, universal bed adapter, wall, pin (16 mm W), (separately available)
ORDER NUMBER	D14197000 - silver



LIGHT DISTRIBUTION

TECHNICAL DRAWING



Photometric data -10% / +20% tolerance

Dcare® **medicool 14 P S3**

TECHNICAL DATA




CARE AND READING LIGHTS

FITTED WITH	1 x compact fluorescent lamp TC-TSE 14W/827 E27 - D68880000
LUMINOUS POWER	450 Lux / 0,5 m
COLOR RENDERING INDEX (CRI)	> 80
COLOUR TEMPERATUR	approx. 2700 K
WORK EQUIPMENT	integrated in spare lamp
CONNECTED LOAD	230-240V; 50Hz
MAINS LEAD	2 m Euro
LUMINAIRE BODY	polycarbonate (PC)
WEIGHT (NET)	approx. 1.5 kg
POWER CONSUMPTION	approx. 14 W
CLASS OF PROTECTION	II
NORMS	EN 60598-1, EN 60598-2-25
SPECIAL FEATURES	separate light strip for orientation
FASTENING	rail, table, universal bed adapter, wall, pin (16 mm W), (separately available)
ORDER NUMBER	D14199000 - silver



LIGHT DISTRIBUTION

TECHNICAL DRAWING



Photometric data -10% / +20% tolerance