


MicroGrip™

Ratcheting Retractor/Grasper for Percutaneous Access

The MicroGrip is a 2.7mm, 28cm long reusable instrument designed to secure and manipulate soft tissue and organs during traditional and Single Port laparoscopic procedures.




Instrumentation & Access in a Single Device

Available in 3 Jaw Patterns

MicroGrip™

Features & Benefits

- Ratcheted Handle
- One Hand Operation
- Nitinol Jaws
- Blunt Dissecting Tip
- Multiple Jaw Configurations
- Eliminates need for additional large trocars and suture closure of incision



Ordering Information:

MicroGrip is available with the following jaw patterns:

Style	Ref:	Description
A	92050-A	Alligator Grasping Jaws
B	92051-B	Bowel Grasping Jaw
C	92052-C	Clutch Grasping Jaw
	92047-RC	Replacement Sheath for MicroGrip, 2.7mm x 28cm
	92049-CB	Cleaning Brush for MicroGrip, 2.3mm dia. X 60cm

For Additional Information Contact Mediflex
or your Local Mediflex Distributor



Mediflex®
SURGICAL PRODUCTS

www.mediflex.com
Toll-free: 1-800-764-0474
Phone: 1-631-582-6424
Fax: 1-631-582-4520