

endo *Look*<sup>®</sup>



- 01 | ARTHROSCOPY
- 02 | GYNECOLOGY
- 03 | LAPAROSCOPY
- 04 | ENT
- 05 | THORACOSCOPY
- 06 | UROLOGY
- 07 | BRONCHOSCOPY
- 08 | PAEDIATRY
- 09 | SCOPE WITH  
ANGLED EYEPIECE**
- 10 | ACCESSORIES

SCOPE WITH ANGLED EYEPIECE

### Scope with Angled Eyepiece

Ø 6,5 mm | 30° | 176 mm

Optik mit Schrägeinblick

Article number <i>Artikelnummer</i>	Technical Data <i>Technische Daten</i>	Description <i>Beschreibung</i>
<b>EN-0650.30.176</b>	Standard Length <i>Standardlänge</i>	176 mm
	Outer Diameter <i>Außendurchmesser</i>	6,5 mm
	Working Channel <i>Arbeitskanal</i>	3,7 mm
	Suctional Channel Diameter <i>Ansaugkanal</i>	1,0 mm
	Irrigation Channel Diameter <i>Spülkanal</i>	1,0 mm
	Direction of View (Degree) <i>Blickrichtung</i>	30°

### Scope with Angled Eyepiece

Ø 5,5 mm | 30° | 200 mm

Optik mit Schrägeinblick

Article number <i>Artikelnummer</i>	Technical Data <i>Technische Daten</i>	Description <i>Beschreibung</i>
<b>EL-0550.30.200</b>	Standard Length <i>Standardlänge</i>	200 mm
	Outer Diameter <i>Außendurchmesser</i>	5,5 mm
	Working Channel <i>Arbeitskanal</i>	1,8 mm
	Suctional Channel Diameter <i>Ansaugkanal</i>	1,0 mm
	Irrigation Channel Diameter <i>Spülkanal</i>	1,0 mm
	Direction of View (Degree) <i>Blickrichtung</i>	30°

