

Sterilization Container



**marSafe**  
*the alternative solution*

## **mar**Safe Sterilization Container Systems *the alternative Solution*

*For more than 3 decades the KLS Martin Group sets standards in hospital sterile supply. The own demand has always been to maximize patient safety from the perspective of packaging systems and to optimize the workflow in the daily practice of CSSD.*

*Thus, we have grown to one of the leading manufacturers of sterilization container systems worldwide. Numerous customers trust in our quality and our experience and therefore choose the sterilization container systems made by KLS Martin.*

### *One brand loaded with unique features*

The marSafe sterilization container system was especially developed to customer needs and desires. The characteristics of the container system were developed based on decades of experience in hospital sterile supply.

### *The new protection cover*

Specifically for the marSafe series a safety cover system was constructed which can be installed or removed at any time. It also protects the sterile barrier from direct exposure. But most important of all: the container can be sterilized with the safety cover. The only permanently flexible protection in the world.

Due to its philosophy “less is more.”

*Less* packing time, *less* storing space, *less* parts to reprocess, **more** efficiency for your CSSD!

### *Reliability? 100%!*

The new closing mechanism guarantees a 100% safe and secure connection between lid and bottom. Two become one – a reliable union.

### *Proven standard for long lasting surfaces*

Unlike other sterilization container systems, we use a so called CompoundPolid surface for the container bottoms. This premium coating makes the surface more resistant towards chemical attack and usual wear and tear.


### *As flexible as you are!*

The variable sealing system – optional side load or top-load – creates great flexibility.

The color coding at the container front allows an easy adaption of colors. Thanks to the new **slide'n go** system! You can eject and insert the color field easily. There are 6 different colors to choose from.

## The KLS Martin quality at a glance

The optional protection cover



Single use paper filter system

Secure stacking, even with the **MicroStop®** family



**slide'n go** – for easy and fast change of the color coding


Seal slot: optional side load or top-load – it's your choice


Bottom with high performance CompoundPolid coating, available with/without perforation


## Sterilization Containers half-size (1/2)




### *marSafe containers, half-size (1/2)*


Front: grey 	Height	closed	perforated
	110 mm	55-230-10-04	55-231-10-04
	140 mm	55-230-13-04	55-231-13-04
	160 mm	55-230-15-04	55-231-15-04
	210 mm	55-230-20-04	55-231-20-04
	270 mm	55-230-30-04	55-231-30-04


Front: red 	Height	closed	perforated
	110 mm	55-240-10-04	55-241-10-04
	140 mm	55-240-13-04	55-241-13-04
	160 mm	55-240-15-04	55-241-15-04
	210 mm	55-240-20-04	55-241-20-04
	270 mm	55-240-30-04	55-241-30-04

Front: blue 	Height	closed	perforated
	110 mm	55-250-10-04	55-251-10-04
	140 mm	55-250-13-04	55-251-13-04
	160 mm	55-250-15-04	55-251-15-04
	210 mm	55-250-20-04	55-251-20-04
	270 mm	55-250-30-04	55-251-30-04


### *marSafe containers, half-size (1/2)*

Front: orange 	Height	closed	perforated
	110 mm	55-260-10-04	55-261-10-04
	140 mm	55-260-13-04	55-261-13-04
	160 mm	55-260-15-04	55-261-15-04
	210 mm	55-260-20-04	55-261-20-04
	270 mm	55-260-30-04	55-261-30-04


Front: green 	Height	closed	perforated
	110 mm	55-270-10-04	55-271-10-04
	140 mm	55-270-13-04	55-271-13-04
	160 mm	55-270-15-04	55-271-15-04
	210 mm	55-270-20-04	55-271-20-04
	270 mm	55-270-30-04	55-271-30-04


Front: light-grey 	Height	closed	perforated
	110 mm	55-280-10-04	55-281-10-04
	140 mm	55-280-13-04	55-281-13-04
	160 mm	55-280-15-04	55-281-15-04
	210 mm	55-280-20-04	55-281-20-04
	270 mm	55-280-30-04	55-281-30-04


## Sterilization Containers three-quarter-size (3/4)




*marSafe containers, three-quarter-size (3/4)*


Front: grey 	Height	closed	perforated
	110 mm	55-232-10-04	55-233-10-04
	140 mm	55-232-13-04	55-233-13-04
	160 mm	55-232-15-04	55-233-15-04
	-	-	-
-	-	-	


Front: red 	Height	closed	perforated
	110 mm	55-242-10-04	55-243-10-04
	140 mm	55-242-13-04	55-243-13-04
	160 mm	55-242-15-04	55-243-15-04
	-	-	-
-	-	-	

Front: blue 	Height	closed	perforated
	110 mm	55-252-10-04	55-253-10-04
	140 mm	55-252-13-04	55-253-13-04
	160 mm	55-252-15-04	55-253-15-04
	-	-	-
-	-	-	

*marSafe containers, three-quarter-size (3/4)*

Front: orange 	Height	closed	perforated
	110 mm	55-262-10-04	55-263-10-04
	140 mm	55-262-13-04	55-263-13-04
	160 mm	55-262-15-04	55-263-15-04
	-	-	-
-	-	-	


Front: green 	Height	closed	perforated
	110 mm	55-272-10-04	55-273-10-04
	140 mm	55-272-13-04	55-273-13-04
	160 mm	55-272-15-04	55-273-15-04
	-	-	-
-	-	-	


Front: light-grey 	Height	closed	perforated
	110 mm	55-282-10-04	55-283-10-04
	140 mm	55-282-13-04	55-283-13-04
	160 mm	55-282-15-04	55-283-15-04
	-	-	-
-	-	-	


## Sterilization Containers full-size (1/1)




### marSafe containers, full-size (1/1)


Front: grey 	Height	closed	perforated
	110 mm	55-234-10-04	55-235-10-04
	140 mm	55-234-13-04	55-235-13-04
	160 mm	55-234-15-04	55-235-15-04
	210 mm	55-234-20-04	55-235-20-04
	270 mm	55-234-30-04	55-235-30-04


Front: red 	Height	closed	perforated
	110 mm	55-244-10-04	55-245-10-04
	140 mm	55-244-13-04	55-245-13-04
	160 mm	55-244-15-04	55-245-15-04
	210 mm	55-244-20-04	55-245-20-04
	270 mm	55-244-30-04	55-245-30-04

Front: blue 	Height	closed	perforated
	110 mm	55-254-10-04	55-255-10-04
	140 mm	55-254-13-04	55-255-13-04
	160 mm	55-254-15-04	55-255-15-04
	210 mm	55-254-20-04	55-255-20-04
	270 mm	55-254-30-04	55-255-30-04

### marSafe containers, full-size (1/1)

Front: orange 	Height	closed	perforated
	110 mm	55-264-10-04	55-265-10-04
	140 mm	55-264-13-04	55-265-13-04
	160 mm	55-264-15-04	55-265-15-04
	210 mm	55-264-20-04	55-265-20-04
	270 mm	55-264-30-04	55-265-30-04

Front: green 	Height	closed	perforated
	110 mm	55-274-10-04	55-275-10-04
	140 mm	55-274-13-04	55-275-13-04
	160 mm	55-274-15-04	55-275-15-04
	210 mm	55-274-20-04	55-275-20-04
	270 mm	55-274-30-04	55-275-30-04

Front: light-grey 	Height	closed	perforated
	110 mm	55-284-10-04	55-285-10-04
	140 mm	55-284-13-04	55-285-13-04
	160 mm	55-284-15-04	55-285-15-04
	210 mm	55-284-20-04	55-285-20-04
	270 mm	55-284-30-04	55-285-30-04

# Sterilization Containers Accessories



## Accessories for marSafe containers

	<b>Item no.</b>	<b>Lid</b>
	55-290-30-04	for container 30 x 30 cm
	55-290-47-04	for container 47 x 30 cm
	55-290-60-04	for container 60 x 30 cm
	<b>Protection cover</b>	
	55-290-99-04	for marSafe container
	<b>Paper filter</b>	
	55-291-01-04	Ø 190/6 mm (500/pack)
	<b>Holder for filter</b>	
	04-250-00-04	for marSafe container
<b>Front</b>		
55-230-00-04	grey	
55-240-00-04	red	
55-250-00-04	blue	
55-260-00-04	orange	
55-270-00-04	green	
55-280-00-04	light-grey	

	<b>Item no.</b>	<b>Seal</b>
	55-996-84-04	for marSafe without indicator (1,000/pack)
	55-996-85-04	for marSafe with indicator (1,000/pack)
	55-996-81-04	for marSafe, topload with indicator (1,000/pack)
	<b>Label</b>	
	55-291-11-04	without indicator, one side adhesive, 35 x 75 mm
	55-291-12-04	without indicator, both sides adhesive, 35 x 75 mm
	55-291-13-04	with indicator, both sides adhesive, 35 x 70 mm
	<b>Coding labels</b>	
	55-291-20-04	no lettering
55-291-21-04	lettered	

## KLS Martin Group

**Karl Leibinger GmbH & Co. KG**  
78570 Mühlheim · Germany  
Tel. +49 7463 838-0  
info@klsmartin.com

**KLS Martin GmbH + Co. KG**  
79224 Umkirch · Germany  
Tel. +49 7665 9802-0  
info@klsmartin.com

**Stuckenbrock Medizintechnik GmbH**  
78532 Tuttlingen · Germany  
Tel. +49 7461 165880  
verwaltung@stuckenbrock.de

**Rudolf Buck GmbH**  
78570 Mühlheim · Germany  
Tel. +49 7463 99516-30  
info@klsmartin.com

**KLS Martin France SARL**  
68000 Colmar · France  
Tel. +33 3 8921 6601  
france@klsmartin.com

**Martin Italia S.r.l.**  
20871 Vimercate (MB) · Italy  
Tel. +39 039 605 6731  
italia@klsmartin.com

**Martin Nederland/Marned B.V.**  
1270 AG Huizen · The Netherlands  
Tel. +31 35 523 4538  
nederland@klsmartin.com

**KLS Martin UK Ltd.**  
Reading RG1 3EU · United Kingdom  
Tel. +44 1189 000 570  
uk@klsmartin.com

**Nippon Martin K.K.**  
Osaka 541-0046 · Japan  
Tel. +81 6 6228 9075  
nippon@klsmartin.com

**KLS Martin L.P.**  
Jacksonville, FL 32246 · USA  
Tel. +1 904 641 7746  
usa@klsmartin.com

**Gebrüder Martin GmbH & Co. KG**  
Representative Office  
121471 Moscow · Russia  
Tel. +7 499 792-76-19  
russia@klsmartin.com

**Gebrüder Martin GmbH & Co. KG**  
Representative Office  
201203 Shanghai · China  
Tel. +86 21 2898 6611  
china@klsmartin.com

**Gebrüder Martin GmbH & Co. KG**  
Representative Office  
Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

**1923**  
**2013**  
**YEARS** | **SURGICAL**  
**INNOVATION**

**Gebrüder Martin GmbH & Co. KG**  
A company of the KLS Martin Group  
Ludwigstaler Str. 132 · 78532 Tuttlingen · Germany  
Postfach 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com

