

***Patient Comfort,
Value, Performance
& Flexibility***

Comfort: Lightweight and compact amplifier comfortable for adults and children

Multipurpose: Support for EEG and PSG video studies in both clinical and ambulatory settings

Standard components: Standard off-the-shelf AA batteries and high speed USB 2.0 connectivity

Performance: Up to 96 hours of ambulatory recording; less than 5 minutes to upload 24 hours of data

Innovation: Unique ambulatory video technology enables up to 48 hours of fully synchronized High-Definition video with simplified setup in patient's home

Efficiency: Record up to 48 hours of High-Definition video fully synchronized over wireless connection

Tools: Microsoft SQL[®] Server database and integrated Spike & Event detection algorithm for maximum efficiency

***Now featuring
High-Definition
Video***





Data Security

Set up individual user or group access and permissions for the system and within the application. Assign separate lock-down rules for different computers on the same network. Access Audit Logs to review user activity on the system.

A High-Performance System with High-Performance Support

Standing behind every XLTEK product is Natus Medical Incorporated, an internationally respected innovator of medical products and services (Nasdaq: BABY).

Our Ambulatory systems are backed by the largest clinical and technical field support network in the Neuro/Sleep industry.

Our in-house support team is staffed with dedicated Technical and Field Service Representatives, as well as dedicated installation Project Managers and Clinical Educators. We provide 24/7 support for LTM, EEG and Sleep Diagnostics, offering the convenience of remote support via WebEx or VPN. Choose from a wide range of extended service contracts to meet your needs, covering software upgrades, preventative maintenance Field Service visits and planned periodic equipment replacement, protecting you against any technological obsolescence.

No one supports you better than Natus.



Natus Medical Incorporated
Excel-Tech Ltd. (XLTEK)
2568 Bristol Circle, Oakville
ON L6H 5S1 CANADA
1-800-303-0306
1-905-829-5300

www.natus.com

P/N 002991B

CE Certified

natus
neurology

Dantec Bio-logic xltek **Stellate**