



Venus-222

22" 24/7 Self-Power Cart Computer



Features

- Intel® J1900 Quad-Core Processor
- Equipped Swappable Batteries x2, support 24hr run-time.
- sBMS Smart Battery Management control tool
- Fully Cable-less, easy deployment
- Aluminum Alloy Rear Cover
- Reading light bar in front
- 22" Full HD high contrast LCD with LED backlight
- Water-Proof : IP65 Front Panel , Rear- IPX1
- 5 Programmable Smart Function Keys
- USB/COM/LAN Medical 4KV Isolation
- High Speed USB 3.0 Ports
- 2.0/5.0MP Pixels Camera (optional)
- 9/12/24V External Power Output (optional)



Specifications

Main Specifications

Processor	Intel® quad-core J1900 (2M Cache, up to 2.42 GHz)
System Memory	DDR3 up to 8GB
Chipset	N/A
OS Support	Microsoft® Windows® 7 (32 bit / 64 bit) Microsoft® Windows® 8.1 (32 bit / 64 bit)
Storage Disk Driver	2.5 inch SATA2 SSD or HDD (optional)
Security	Trusted Platform Module, Smart Card Reader, RFID (optional)
Wireless Communication	Intel 802.11 a/b/g/n (optional) Bluetooth (optional)
Speaker	2W x 2
Camera	2.0/5.0MP Pixel (optional)
Function Key	Power On/Off, LCD Brightness Up, LCD Brightness Down, Touch Screen On/Off, Reading Light On/Off
Remote Control	IR Receiver (optional)

Display

Size	22 inch VA
Resolution	1920 x 1080
Max. Colors	16.7M
View Angle	178°(H) / 178°(V)
Contrast Ratio	3000:1
Brightness	250nits
Touch Screen	5-Wire Resistive Touch Screen
Sensor	Light Sensor

I/O

USB	4KV Isolated USB 2.0 x 1 , USB 2.0 x 4, USB 3.0 x 2
COM	4KV Isolated RS-232 x2
Ethernet	4KV Isolated Gigabit LAN x 1
Video Out	Display Port x 1
Audio	Mic-in x 1, Audio-out x 1
External Power Output	9/12/24V rear power output (optional)

Battery Pack

Battery Type	6 Cell Lithium Battery
Battery Capacity	5800 mAh
Battery Output Voltage	10.8V
Battery Charge Time	2 hrs for 85%, 3hrs for 100%
Battery Run time	6hrs *For more details, please see the user guide

Mechanical and Environmental

Operating Temperature	0°C ~ 40°C(32°F ~ 104°F)
Storage Temperature	-20°C ~ 60°C(-4°F ~ 140°F)
Storage Humidity	10%~95% @ 35°C, non-condensing
Degree of Protection	Front Panel IP-65 Certified, IP x 1
System Cooling	Fan-less
EMC	CE/FCC Class B, UL-60601-1, EN-60601-1,PSE, RCM