


Unlimited possibilities
and quick positioning of
the patient by four power-
ful motors. Multifunctional
arm-, head-, and foot-rests
meet all wishes. Conveni-
ent: included castor faci-
lity for effortless movement.



PROMAT MX. Diversity by details.



MEDICAL

PROMAT MX. The technical data.

Part	Max. load	Adjusting				Dimensions (lengthxwidth)
		Range	Angle	Speed unloaded	Speed loaded	
Head rest	70 kg	100 mm	180°	manual	manual	240 x 375 mm
Back rest		-	0° - 78°	14 sec	20 sec	590 x 545 mm
Seat	70 kg	-	0° - 20°	9 sec	9 sec	500 x 540 mm
Leg rest	35 kg	-	0° - 78°	9 sec	12 sec	310 x 530 mm
Foot rest		160 mm	180°	manual	manual	195 x 445 mm
Height adjustment	175 kg	310 mm	-	14 sec	22 sec	from 595 - 905 mm
Total weight	approx. 65 kg					
Total length	1820 - 2160 mm					
Total width	545 mm (lying surface) / 745 mm (arm rests)					



Foot rest adjustment



Foot control



Multifunctional remote control

PROMAT MX



PROMAT MX. The accessories.

- Wide armrests
- Multifunctional headrest
- Music (CD/Speaker)
- Lights (under the bottom frame covering)
- Seat and back rest heating
- Various cover materials
- Foot control



Wide armrests



Multifunctional headrest