

AEROCASE® - Pro 1R First Responder

Emergency bags | Emergency backpacks



new
HUM - NOVELTY

First Responder

Reliable companions for the comprehensive first aid

AEROcase[®] - Pro

1R First Responder

Others speak of quality, we have quality!

AEROcase[®] - Pro^{1R} emergency bags and backpacks provide the safe transport of equipment and are reliable partners when it comes to first aid.

By the compact and lightweight design they are easy and convenient to transport and ideal companion for locations with long distance or restricted room.

Easy accessible outside pockets and (useful interiors for) well structured interiors allow quick access to the equipment on site. Highly fluorescent reflective strips provides for safety during day and night.


Due to the high-quality materials and manufacturing, the emergency bags and backpacks are extremely durable and well equipped for the high loads in everyday use.

Possible areas of application

- Medical service
- Ambulance
- Special forces
- Civil protection
- Companies
- Schools
- Private

AEROcase[®] - Pro^{1R} Material details


Materials for highest requirements and durability



AEROTex[®] - Dura

AEROTex[®] - Dura stands for high quality polyester (600 D).

The nylon fibers are PVC-coated on the inner surface and are extremely robust and give the material the necessary properties such as tear- and abrasion-resistance for daily use.



AEROTex[®] - Plan

AEROTex[®] - Plan is a high quality and robust plane material.


It was especially designed for professional use in emergency medical service. Due to the both-sided PVC-coating, the plane material is very resistant to loads and is furthermore waterproof. It is characterized by scratch-, tear- and abrasion-resistance and is in addition also very hard-wearing. The surface can be cleaned with appropriate cleaners and also with appropriate disinfectants.

AEROcase® - Pro^{1R} Processing details

The common details at one glance

Handles

The ergonomically formed plastic handles with reinforced core and the soft grip handles with hook and loop fastener provide a high comfort. The eye-catching colouring also ensures an improved cognition and a fast access when time is short.



Zipper

All zippers are provided with coloured zipper extensions. The opening of the pockets and compartments can be simplified and thanks to the improved cognition a fast access to the equipment is possible.



Transparent pocket

An additional clear window for a good visible placement of documents and information, is placed on the outside.




Hook and loop strip for patches

A hook and loop strip is fixed on the outside to attach individual patches. (Exception: AEROcase® - PRO^{1R} BS1)




Highly fluorescent reflective strips

The large reflective stripes (yellow, silver gray) with a high graded light reflection ensure optimum visibility and safety on site.



Bottom area made of plan material

The bottom area is made of black, durable and washable plan material. It is raised trough shaped, so that the lower edges of the outer surfaces are also coated. The advantages are: Protection against humidity, dirt and damage as well as easy cleaning.



Pyramid feet

The pyramid feet are made of very compressible plastic and give the emergency bag and backpack high stability and also protect the bottom area against abrasion.

* Only for models made of black polyester

AEROcase® - Pro1R PL1C



Versatile applications, high functionality and high-quality processing qualify the Pro1R PL1C for the perfect standard emergency backpack for professional use. The interior impresses with its flexibility and individual organizational possibilities. In the three divided main compartment, which can be closed by a cover flap, the integrated 3 matching module pockets can be integrated easily. The additional 2 modular pockets can be placed on the front with hook and loop fastener. Additionally 2 side pockets and 2 front pockets provide storage space for extensive medical equipment. The 5 reflective strips ensure an all around visibility. Extra wide and padded shoulder straps with stomach and upper body belts ensure a comfortable fit of the backpack.

Polyester

36 x 47 x 26 cm | 2 | ca. 23

Plane

36 x 47 x 26 cm | 3 | ca. 23



Ergonomic formed handles

- Grip comfort due to the reinforced core
- Better perception of signal colours

Transparent pocket

- For good visible placement of documents and information

2 Front pockets

- Upper one for writing material
- The one below with plenty of space for ampoule cases | intubation cases etc.
- Additional zipper compartment for documents

Wide hook and loop fastener

- To fix individual patches

5 Reflective strips

- Front pocket | side pockets | carrying strap
- All around visible and safe

2 Side pockets

- Space for blood pressure meter / stethoscope / sharps container etc.

Conspicuous zipper

- Better cognition
- Easy and fast opening

AEROCASE® - Pro1R PL1C

Features and models



Clear arranged interior concept

The modular and structured interior concept ensures a secure placement of the equipment.

- Fold-out front
- Main pocket contains removable frame with extra separating strip
- 5 module pockets for interior:
 - 3 in corresponding divided main compartment
 - 2 can be placed inside the front pocket
- Elastic attachment straps (quality A)
- Lockable main pocket with cover flap with hook and loop fastener

Protected all around

- Dimensionally stable exterior walls
- Protected bottom area

Modell: Plane

Modell: Polyester



Individual organization

- 5 module pockets in different sizes
 - Plane: plain-coloured (red or blue)
 - Polyester: assorted colours
- transparent and with handles

Comfortable carrying system

- Padded shoulder straps
- Additionally quick release belts for fixation on stomach and upper body

Polyester | AEROTex® - Dura

Plan | AEROTex® - Plan



HT01-RPL1C-R



HT01-RPL1C-B



HT01-RPL1C-BLK



HT03-RPL1C-R



HT03-RPL1C-B

AEROCASE® - Pro1R BS1



is the smallest model among the Pro1R emergency bags and it is ideal for storing small basic equipment. The storage space is divided into the actual main compartment with a divided interior as well as 1 large front pocket and 2 large side pockets. Security on site thanks to the 3 highly fluorescent reflective strips and a comfortable fit are standard. Useful details, such as the padded and removable inner partition, the area which is divided into 2 sectioned front pockets, the comfortable carrying handle as well as the system with quick-release buckle and hook and loop fastener makes the "small" one to an even "bigger" one.

Polyester

37 x 21 x 27 cm | 0,7 kg | ca. 12

Ergonomic formed handles

- Grip comfort due to the reinforced core
- Better cognition because of signal colours

Large transparent pocket on the backside

- Ideal for the good visibility of documents and information

Shoulder carrying strap

- Adjustable in length
- High quality spring clasp

Quick-release buckle and hook and loop system

- Quick release belt on the front
- Hook and loop fastener on the front cover
- Fast and easy access to the interior



2 Large side pockets

- Lots of optional storage room

Sectioned front pocket

- Big front pocket with divider
- More sorting options

Conspicuous zipper

- Better cognition
- Easy and fast opening

3 Reflective strips

- Front pocket | side pockets
- All around visible and safe

AEROCASE® - Pro^{1R} BS1

Features and models



Padded interior layout

- Protects the equipment
- Ensures order



Flexible

- Fully removable interior layout
- Individual

Clear interior concept

- Zipper compartment within the cover
- Elastic attachment straps on the cover as well as on the front
- Divided front pocket below
- Padded compartments, including pockets
- Side clips protect the interior from dirt and moisture



More room for maneuver

- Delivery including shoulder strap
- Free hands for more mobility

Transparent pocket on the entire back

- Much space to place documents and information for easy perceptibility

Polyester | AEROTEX® - Dura



HT01-RBS1-R



HT01-RBS1-B

AEROcase® - Pro^{1R} BM1

1R First Responder

The compact Pro^{1R} BM1 is made for sophisticated users and has additional space for expanded basic equipment. The main pocket with its three-part, related to practice interior is also supplemented with 2 spacious side pockets. It is easy to open and allows a quick access to the clearly structured interior within seconds. The construction and durable material ensures dimensional stability and protection against mechanical damage. The 4 highly fluorescent reflective strips on the exterior sides provide optimal safety. The 4 highly fluorescent reflective strips on the exterior sides provide optimal safety. The padded hook and loop handle and the included shoulder strap ensure wearing comfort.

Polyester

43 x 25 x 27 cm | 1 kg | ca. 18

Plane

43 x 25 x 27 cm | 2 kg | ca. 18

Large transparent pocket on the cover

- Ideal for the good visibility of documents and information

Hook and loop handle

- Soft padded handle
- For comfortable carrying

Conspicuous zipper

- Better cognition
- Easy and fast opening



2 Large side pockets

- Generous storage room
- With robust zippers

4 Reflective strips

- Side pockets | front side | back side
- Very good light reflection
- All around visible and safe

Protected bottom area

- Trough shaped bottom made of stable plan material
- 4 robust plastic pyramid feet
- Stability and protection

AEROCASE® - Pro1R BM1

Features and models

Flexible

- Kompaktes und platzsparendes Design
- Multifunctional storage compartments
- Easy to use



Large transparent pocket

- Much space for documents and information to identify quickly



Protected all around

- Very robust and durable
- Dimensionally stable construction
- Protected bottom area



More room for maneuver

- Delivery including shoulder strap
- Free hands for more mobility

Interior concept

- Large and clear space for extended basic equipment
- Zipper compartment within the cover
- Main compartment divided into three sections with fixed internal division
- Reasonable arrangement of elastic attachment straps (quality A):
In the cover | At the inner back side



Polyester | AEROTex® - Dura

Plan | AEROTex® - Plan



HT01-RBM1-R



HT01-RBM1-B



HT03-RBM1-R



HT03-RBM1-B

AEROCASE® - Pro^{1R} BL1



Die AEROCASE® Pro^{1R} BL1 is the ideal emergency bag for professional first aid. It is extremely robust, durable and dimensionally stable constructed and equipped with wide and high quality zippers. The padded walls protect the entire equipment. The space concept is dimensioned for the specific requirements of an extensive and well-organized first aid. The bag is featured with a large main compartment with a variable interior and numerous large pockets at the sides, the front and on the cover. The hinged front pocket with its many useful elastic attachment straps and compartments is divided extensively. The shoulder strap and the additional back strap in a pocket on the back ensure a versatile wearing comfort.

Polyester

52 x 32 x 32 cm | 1,6 | ca. 28

Plan

52 x 32 x 32 cm | 2,7 | ca. 28

Large cover pocket

- Additional, optimal storage room

Transparent pocket

- For good visible placement of documents and information

5 Reflective straps

- Front side | back side | side pockets | cover
- A very good light reflection
- All around visible and safe



2 Side pockets

- With robust zippers
- With additionally 2 large elastic attachment straps (right)

Conspicuous zipper

- Better perception
- Easy and fast opening

Large front pocket

- With many useful compartments and elastic attachment straps

Robust bottom area

- Trough shaped bottom made of stable plan material
- 4 robust plastic pyramid feet
- Stability and protection

AEROCASE® - Pro1R BL 1


Features and models

Professionally organized

- Very large space for extensive equipment
- Zipper compartment within the cover
- Structured, practical interior layout
- Interior with variable divider system
- Robust elastic attachment straps (quality A) in:
 - Interior
 - Side pocket
 - Front pocket

Variable interior

- Removable
- Movable



Shoulder strap

Optimum protection

- Reinforced outer walls
- Completely padded
- Protected bottom area

Divided front pocket

- Fold-out front
- Many useful straps of different sizes and dividers
- Useful e.g. for:
 - Infusion accessories
 - Diagnostic equipment
 - Etc.



Back straps

Versatile wearing options

- Back straps in the pocket on the backside
- Additional shoulder strap
- Free hands for more mobility

Polyester | AEROTex® - Dura

Plan | AEROTex® - Plan



HT01-RBL1-R



HT01-RBL1-B



HT01-RBL1-BLK



HT03-RBL1-R



HT03-RBL1-B

AERObase® - Pro^{1R} PL1C

Item No.	Description	Size	Material	Colour	Dimensions *	Weight	Volume
HT01-RPL1C-R	AERObase® - Pro ^{1R} PL1C Emergency Backpack	L	Polyester	Red	36 x 47 x 26 cm	2 kg	ca. 23 L
HT01-RBPL1C-B	AERObase® - Pro ^{1R} PL1C Emergency Backpack	L	Polyester	Blue	36 x 47 x 26 cm	2 kg	ca. 23 L
HT01-RPL1C-BLK	AERObase® - Pro ^{1R} PL1C Emergency Backpack	L	Polyester	Black	36 x 47 x 26 cm	2 kg	ca. 23 L
HT03-RPL1C-R	AERObase® - Pro ^{1R} PL1C Emergency Backpack	L	Plan	Red	36 x 47 x 26 cm	3 kg	ca. 23 L
HT03-RBPL1C-B	AERObase® - Pro ^{1R} PL1C Emergency Backpack	L	Plan	Blue	36 x 47 x 26 cm	3 kg	ca. 23 L

*dimensions incl. front and side pockets



AERObase® - Pro^{1R} BS1

Item No.	Description	Size	Material	Colour	Dimensions *	Weight	Volume
HT01-RBS1-R	AERObase® - Pro ^{1R} BS1 Emergency Bag	S	Polyester	Red	37 x 21 x 27 cm	0,7 kg	ca. 12 L
HT01-RBS1-B	AERObase® - Pro ^{1R} BS1 Emergency Bag	S	Polyester	Blue	37 x 21 x 27 cm	0,7 kg	ca. 12 L

*dimensions incl. front and side pockets



AERObase® - Pro^{1R} BM1

Item.No.	Description	Size	Material	Colour	Dimensions *	Weight	Volume
HT01-RBM1-R	AERObase® - Pro ^{1R} BM1 Emergency Bag	M	Polyester	Red	43 x 25 x 27 cm	1 kg	ca. 18 L
HT01-RBM1-B	AERObase® - Pro ^{1R} BM1 Emergency Bag	M	Polyester	Blue	43 x 25 x 27 cm	1 kg	ca. 18 L
HT03-RBM1-R	AERObase® - Pro ^{1R} BM1 Emergency Bag	M	Plan	Red	43 x 25 x 27 cm	2 kg	ca. 18 L
HT03-RBM1-B	AERObase® - Pro ^{1R} BM1 Emergency Bag	M	Plan	Blue	43 x 25 x 27 cm	2 kg	ca. 18 L

*dimensions incl. front and side pockets



AERObase® - Pro^{1R} BL1

Item.No.	Description	Size	Material	Colour	Dimensions *	Weight	Volume
HT01-RBL1-R	AERObase® - Pro ^{1R} BL1 Emergency Bag	L	Polyester	Red	52 x 32 x 32 cm	1,6 kg	ca. 28 L
HT01-RBL1-B	AERObase® - Pro ^{1R} BL1 Emergency Bag	L	Polyester	Blue	52 x 32 x 32 cm	1,6 kg	ca. 28 L
HT01-RBL1-BLK	AERObase® - Pro ^{1R} BL1 Emergency Bag	L	Polyester	Black	52 x 32 x 32 cm	1,6 kg	ca. 28 L
HT03-RBL1-R	AERObase® - Pro ^{1R} BL1 Emergency Bag	L	Plan	Red	52 x 32 x 32 cm	2,7 kg	ca. 28 L
HT03-RBL1-B	AERObase® - Pro ^{1R} BL1 Emergency Bag	L	Plan	Blue	52 x 32 x 32 cm	2,7 kg	ca. 28 L


*dimensions incl. front and side pockets



AERObase® - Pro^{1R} Innentaschen-Sets

Item.No.	Description	Material	Colour
HT01-RPL1-C-SET	AERObase® - Pro ^{1R} Inside Pocket Set	Polyester	5 farbig
HT03-RPL1-C-SET-R	AERObase® - Pro ^{1R} Inside Pocket Set	Plan	Red
HT03-RPL1-C-SET-B	AERObase® - Pro ^{1R} Inside Pocket Set	Plan	Blue

*dimensions incl. front and side pockets



HUM

HUM Gesellschaft für Homecare
und Medizintechnik mbH
Zum Pier 79 | D-44536 Lünen

FON +49 (0)2 31/ 88 08 85-0 | FAX +49 (0)2 31/ 88 08 85-58
Mail: sales@hum-online.de | http://www.hum-online.de

AERObase® - Pro^{1R}
PP-AERObase_Pro^{1R}_GB | Stand 16.06.2014 | Revision 2
© HUM Gesellschaft für Homecare und Medizintechnik mbH