

Portable Surgery Diode Laser Systems



Wuhan Gigaa Optonics Technology Co.,Ltd.

www.gigaalaser.com

Portable Surgery Diode Laser Systems

SOFT TISSUE SURGERY LASER

The new Gbox portable surgery diode laser system from GIGAA is the best value on the market today. Small size body case with 15W power, this mini surgery laser system has 810nm and 980nm dual-wavelengths, 810nm/980nm with pilot beam 300mW for choice, it can be used in the continuous as well as pulsed mode with a contact or non-contact handpiece.

Features

- Compact design makes it small in size and light in weight, easy to transport;
- Big color touch screen, Encoder knob dual operation;
- Highly integrated and modular design;
- New design inside to make maintenance easy and low cost;
- Advanced software with presetting protocols to be manually customized by operator;
- Built-in power meter can test the power in the software by user;
- Output power calibration at the end of fiber to make the maintenance easy and low cost;
- Adopt air cooling system guarantee stable performance and long life time.



Main Clinical Applications

- Dentistry
- Lipolysis (Laser fat burning)
- EVLT (Endovenous laser treatment)
- Spider Vein
- ENT
- Dermatology
- Orthopaedics
- PLDD (Percutaneous Laser Disc Decompression)
- Neurosurgery
- Gynaecology
- General surgery
- Therapy
- vet



Clinic cases

Please visit our website www.gigaalaser.com for more videos.

Accessories

Standard



Dental Handpiece




Fiber



Fiber stripper



Safety Goggle



Fiber Cutter



Foot Switch


Optional



Therapy Handpiece



Bleaching Handpiece



ENT Handpiece



Lipolysis Handpiece



Focus Handpiece



ADD: B8-A5, Hi-Tech Medical Device Industrial Park, #818 Gaoxin Avenue,
Wuhan 430206, China
Ph: 86-27-67848871 67848872 Fax: 86-27-67848873
<http://www.gigaalaser.com> Email: info@gigaalaser.com



CHEESE™ Velas™



Portable Surgery Diode Laser Systems



532nm Green Laser

Surgery
DEEP Tissue Therapy
Veterinary
Beauty (3W/5W/8W-532nm Green Laser)

Technical Specifications

GBOX

Type	GBOX-10B/C/F	GBOX-15A/B	GBOX-15AB	GBOX-15AP/BP	GBOX-3-8G/GC
Wavelength	980nm/940nm/1064nm	810nm/ 980nm	810nm+980nm	810nm/ 980nm+635nm	532nm/532nm+940nm
Maximum Power	10W	15W	15W+15W	15W+300mW	3W-8W/5W+10W
Operation Mode	CW, Single Or Repeat Pulse				
Pulse Duration	25µs-10s				
Repetition Rate	0.05Hz-20KHz				
Pilot Beam	Red Diode Laser Of 635nm, Power < 5mW				
Control Mode	Touch Screen, 8" True Color				
Dimension	315(H)*215(W)*245(L)mm				
Lifetime	Over 10000 hrs				
Warranty	12 months free warranty, 5 years technical support				
Weight	4Kg				
Packing	Carry case with carton				