



ECOspin III

Compact Table Top Centrifuge
for daily routine

12-fold fixed-angle rotor for 15ml tubes
Brushless Motor with electronic speed
control by 100 r/min

Digital adjustment of operating time and rotating speed
Hygienic membrane keyboard
Electro-magnetic lid lock
Tough, easy-to-clean plastic case



medical **ECONET**

GERMANY

Spezifikationen ECOspin III

Equipped with a 12 - angle rotor for tubes with a volume of 15 ml and a max. RCF of 2.400 g, ECOspin III is the optimal device for a general practitioner.

The ECOspin III centrifuge needs just little storage place and is very easy to use. The last selected setting for rotational speed and time can be saved. It is ready to use at the touch of the button.

For users safety, the lid gets locked electro-magnetically during operation. The durable plastic case and hygienic membrane keyboard enable an easy cleaning.



Dimensions

255 mm (W) x 245 mm (D) x 140 mm (H)



Performance data

Rotor	12 x 15ml
Rotation frequency	increase by 100 min ⁻¹
Max. RCF (relative centrifugal force)	2.400 g
Max. Speed	4500 r / min
Min. Speed	300 r / min
Power	110-240V 50-60 Hz
Noise	ca. 56 dB (A)
Dimensions	255 x 245 x 140mm
Weight	6 kg

medical **CONET**
GERMANY

Im Erlengrund 20, 46149 Oberhausen - Germany

t: +49 (0)208 / 377 890 - 0 info@medical-econet.com
f: +49 (0)208 / 377 890 - 55 www.medical-econet.com

distributor

*All specifications are subject to change without notice.