



ECOtom 80

ECOtom 80 / 200 / 400

Efficiency and precision

Stable and powerful generators
Highest patient safety through REM system
Digital value indicator of chosen power output
Monopolar and Bipolar applications
Activation through Hand- or Footswitch
Safe application through color-coded membrane keypad



medical **ECONET**

GERMANY

Specifications

ECOtom	80	400	200
Max. Output	80W	400W	200W
Pure cut	80/300Ω	400/300Ω	200/300Ω
Blend	80/300Ω	250/300Ω	150/300Ω
Blend II		200/300Ω	120/300Ω
Blend III		150/300Ω	100/300Ω
Contact Coagulation	50/300Ω	120/300Ω	100/300Ω
Spray Coagulation		100/300Ω	80/300Ω
Bipolar Coagulation	70/300Ω	99/100Ω	80/300Ω
Power Frequency	1.6 MHz	400 KHz	
Power source	115 - 230V, 50 - 60 Hz	115 - 230V, 50 - 60 Hz	
Safety class	I BF	I CF	
Neutral electrode	- NF -	- F -	
Weight	6 Kg	13 Kg	
Dimensions	270x310x110 mm	330x455x150 mm	

ECOtom - Series

Power stable generators for monopolar and bipolar applications make ECOtom to an ideal appliance for operating rooms, outpatient clinics and practitioners.


HF - efficiency can be finely set up by 1W increments, allowing a precise adjustment of HF - power to the requirements during surgery

Accessories


■ Bipolar Forcep



■ Disposable Double-Handswitch




■ Monopolar-Handswitch




■ Neutral (silicon) Electrode



■ Patient Return Plate Cable



■ Ground Cable



■ Bipolar Footswitch



■ Double Footswitch



medical **ECONET**
GERMANY

Im Erlengrund 20, 46149 Oberhausen - Germany

t: +49 (0)208 / 377 890 - 0 info@medical-econet.com
f: +49 (0)208 / 377 890 - 55 www.medical-econet.com

distributor

*All specifications are subject to change without notice.