The ECG with more knowledge built in

Introducing the MAC* 5500 HD





If you can't see the pacemaker, you're not seeing the whole patient

During a patient's cardiac event, physicians have moments to make a diagnosis. Modern pacemakers are difficult to detect in standard ECGs, and inferring their presence takes time many clinicians simply don't have in an emergency situation.

Missed or inaccurate pacemaker detection can lead to incorrect diagnosis, wasted treatment time and costs, and increased risk for patients.

The MAC 5500 HD incorporates new tools for ECG acquisition, enabling confident clinical decisions with a three-step solution that delivers on speed and accuracy.



Advancing clinical excellence in three steps

The MAC 5500 HD combines proven 12SL* analysis with technology advancements to provide one of the most advanced ECG pacemaker detection systems available.

- **Detect pacemaker pulses more accurately**¹ via a high sample rate, reducing the risk of treatment contraindicated for paced patients.
- Advanced interpretation of paced rhythms via the Marquette* 12SL program. The program can detect and report the underlying rhythm, increasing the breadth of decision support capabilities.
- Clear visualization and annotation of pacer pulses through MUSE* v8 Cardiology Information System, which provides a dedicated pace-annotation channel. This helps to reduce the need to identify pacer pulses within the ECG signal, supporting efficient diagnosis.


The latest edition of Marquette 12SL provides analysis for detecting bi-ventricular pacemakers, identifying

the underlying rhythm, as well as the chamber(s) being paced.

| Internal Content of the Content

Proven power: HD Pacemaker workflow

Clear data is better data. The CAM HD acquisition module reduces noise and artifact for clearer ECG tracings. This enhanced clarity improves the system's ability to accurately detect the presence of pacemaker pulses and can decrease the need to redo ECGs due to poor data quality.



When seconds count, workflow matters. The HD Pacemaker feature fosters a seamless flow of information from the MAC 5500 HD and into MUSE v8.







New capabilities, familiar features, well known user interface

The MAC 5500 HD system brings you improved functionality and clinical workflow while maintaining the familiar look and ease of use you expect from a GE Healthcare ECG system.

• Familiar look and feel - GE's familiar user interfaces between the different ECG devices can help minimize the training needs for your staff.



©2011 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

*MAC, Marquette, and MUSE are registered trademarks of General Electric Company.

*12SL and MobileLink are trademarks of General Electric Company.

GE Healthcare, a division of General Electric Company

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement, and performance solutions services helps our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our "healthymagination" vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access, and improving quality around the world. Headquartered in the United Kingdom, GE Healthcare is a unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employees are committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com.

GE Healthcare P.O. Box 900, FIN-00031 GE, Finland Tel. +358 10 394 11 Fax +358 9 146 3310

www.gehealthcare.com

