

BELLON *energy efficient and ecologic*

Examination lamp or for minor surgery intended for cares units and private doctors wishing a powerful and precise lighting.

- 17 W
- 5500° Kelvins
- 1440 Lumens
- Guarantee 3 years
- LED's lifetime: 50 000 hours
- Color Rendering Index (IRC) 95
- Positioning of the articulated arm compensated with spring, incorporated into the joint (soft and precise)
- Orientation and positioning of the head facilitated by the flexible ending (L 140 mm).
- No contact switch on the trolley, actionable by foot without any risk of inadvertent activation.
- Double intensity switch on the head of the lamp
- High capacity trolley basis, 4 brake-equipped casters, very stable and very handy.
- Footprint : 720 x 460 mm.
- Vertical Matt L 1130 mm, height adjustable (+/-) 300 mm.
- Directional arm L 900 mm.
- Total weight : 20 kg.
- Reinforced carton packaging: 21 x 60 x 120 mm.
- Installation time : 3 - 4 minutes.
- Cable length: 4 m
- Insulation class I
- Medical device class I
- at 30 cm: 55000 lux, illuminated field: 8 cm
- at 50 cm: 42000 lux, illuminated field: 10 cm
- at 100 cm: 13000 lux, illuminated field: 20 cm
- Code: LED17900DI-ISC

