

DIANA Halogen lamp

Powerful and versatil lamp, with an ergonomic handle



- 50 watts
- 3000° kelvins
- 1000 Lumens
- Guarantee 3 years
- Bulb's lifetime 5000 hours (Dichroic source GU 5,3 OSRAM Decostar ES SP 10° high performance)
- Color Rendering Index (IRC) 100
- Cable length 3 meters allowing maximum mobility
- Insulation Class II
- Available colors: white, black
- Lamp's head dimensions 110 x 110 mm
- Lamp's base dimensions 69 x 69 mm, H: 115 mm, inter axial: 54,5 mm
- Distance 30 cm: 55000 lux, Illuminated field: 5 cm
- Distance 50 cm: 31000 lux, Illuminated field: 7 cm
- Distance 100 cm: 15000 lux, Illuminated field: 17cm
- Flexible length: 650 mm or 1000 mm
- Switch on/off, dimmer or no contact switch
- Option with battery on the base: 2 h autonomy (charger supplied)
- Ease of use 2/3 less heat.
- Code: BOR50651 – BOR51001 – BOR50651V – BOR51001V
BOR50651-ISC – BOR51001-ISC – BORBAT50651V
BORBAT51001V



lid-lamps.com