

HEPTA Neutral white LED

Energy efficient and ecologic 

Very polyvalent examination lamp, powerful and comfortable, without any heat output, optimum handling



- 7 watts
- 4000° kelvins
- 900 Lumens
- Guarantee 3 years
- LED lifetime 50 000 hours (Cold neutral white light, no UV or infrared output)
- Color Rendering Index (IRC) 95
- Cable length 3 meters allowing maximum mobility
- Insulation Class II
- Available colors: white, black
- Lamp's head dimensions 110 x 110 mm
- Lamp's base dimensions 69 x 69 mm, H: 115 mm, inter axial: 54,5 mm
- Distance 30 cm: 53000 lux, illuminated field: 10 cm
- Distance 50 cm: 30000 lux, illuminated field: 15 cm
- Distance 80 cm: 14000 lux, illuminated field: 22 cm
- Distance 100 cm: 9500 lux, illuminated field: 25 cm
- Flexible length: 650 mm or 1000 mm
- Switch on/off, dimmer or no contact switch
- Option with battery: 2h autonomy (charger supplied)
- Code : LED07651-LED07651V-LED071001-LED071001V
LED06BAT650-LED06BAT1000

CE



Company certified ISO 9001 by



lid-lamps.com