

# LITA Halogen reading lamp

Reading lamp that fits over the gallows of medical beds through its clamp.



- 230/12V/20 watts. Halogen bulb with large beam (36°) for a comfortable reading
- 3000° Kelvins
- 400 Lumens
- Guarantee 3 years
- Bulb's lifetime 5000 hours (Dichroic source GU 5,3 OSRAM)
- Color Rendering Index (IRC) 100
- Flexible and omnidirectional arm: 250 mm
- Cable length 3 meters allowing maximum mobility
- Pull switch on the head of the lamp for easy use for bedridden person.  
(+ Second switch on the base of the lamp)
- Insulation Class II. Built-in transformer in the base lamp
- Available colors: white, black
- Lamp's head dimensions 110 x 110 mm
- Lamp's base dimensions 69 x 69 mm, H: 115 mm, inter axial: 54,5 mm
- Ease of use 2/3 less heat.
- Code : LIT20250+ETRPM006



Company certified ISO 9001 by



lid-lamps.com