

THERA Infrared lamp

Heating lamp, ideal for physiotherapy, massages, osteopathy, acupuncture...



- 250 W
- Wavelength 1200 nm
- Guarantee 3 years
- Bulb's lifetime IR (E27) : 5000 hours (integrated reflector in hard glass bulb, 90% of the energy is transmitted as an infrared)
- Cable length 4 meters allowing maximum mobility
- Insulation Class I
- Available colors: white
- Lamp's head diameter 260 mm
- Lamp's base dimensions 69 x 69 mm, H: 115 mm, inter axial: 54,5 mm
- Flexible length: 860 mm
- Switch on/off or no contact switch
- Improves recovery after exercise, by a natural and safe method.
- Recommended for specialties such as: physiotherapy/massage/acupuncture/orthopedics.
- Robust construction, easy handling, accurate and reliable positioning
- Protection grill with three silicone covered fixations points.
- Ergonomic handle
- Code : IRT250860 – IRT250860-ISC



lid-lamps.com