

VERA *energy efficient and ecologic*



LED magnifying lamp, very appreciated in dermatology and for all precise examinations



- 4 W
- 6000° Kelvins
- 300 Lumens
- Guarantee 3 years
- LEDs' lifetime 50000 hours
- Color Rendering Index (IRC) 95
- Cable length 1,5 meter
- Insulation Class II
- Color: white
- Diopter: 5
- Magnifying glass diameter 120 mm, magnification: x 2.25
- Length of the articulated arm: 860 mm
- Clamp provided
- Made in PRC
- Code: LED00912
- Mounting accessories : heavy trolley basis PIR00912, wall bracket EMU00912, universal clamp ETRPM006, clamp for medical rails ERTPM007



CE

ISO 9001
BUREAU VERITAS
Certification



lid-lamps.com