

# Flowstar

## Uroflowmeter



*The Flowstar is a portable digital flowmeter designed for practical, everyday flow study investigations, where ease of use, simplicity and cost effectiveness are top priorities.*

### ■ Automatic Investigation and Analysis

Investigation automatically starts when patient begins voiding and directly prints flow results.

### ■ Digital Flowsensor

Weighted transducer type flowsensor, integrated in an adjustable stand.

### ■ Customizable setting

Flow and volume scales can easily be adjusted. Pre- and post- flow times are user selectable.

### ■ Siroky Flow Plots

Siroky flow plots can be printed to diagnose normal or obstructed urinary flow in male patients.

### ■ Compact

Flowstar is small, lightweight and portable. System can easily be set up and used in any environment.



Medical Measurement Systems

P.O. Box 580

7500 AN Enschede

The Netherlands

T + 31 - 53 - 480 37 00

F + 31 - 53 - 480 37 01

E [info@mmsinternational.com](mailto:info@mmsinternational.com)

I [www.mmsinternational.com](http://www.mmsinternational.com)