

Orion II

Ambulatory pH Recorder



The Orion II is without doubt one of the most comprehensive user friendly pH recorders available to diagnose in Gastro-Esophageal Reflux Diseases (GERD) and non cardiac chest pain.

■ Recording up to 4 pH channels

Glass and Antimony with LES locator.

■ Virtual Instructor Program™ (VIP), very easy to use

Guides user to successfully complete procedures with minimal training.
Easy calibration of catheter via the PC.

■ On-line View of Recordings

Wireless Bluetooth® techniques make it possible to view data while recording.
High quality signal recording can be guaranteed before patient leaves investigation room.

■ Small and Lightweight

Compact housing with comfortable carrying bag. Extremely small and lightweight for maximum patient comfort. A pediatric bag is also available.

■ Sleeplab (optional)

A sleeplab connector can be added as an option, for specs ask MMS



Medical Measurement Systems

P.O. Box 580

7500 AN Enschede

The Netherlands

T + 31 - 53 - 480 37 00

F + 31 - 53 - 480 37 01

E info@mmsinternational.com

I www.mmsinternational.com