

Orion II SL

Ambulatory pH recorder with Sleeplab connection




The Orion II SL can easily be mounted to a Sleeplab System. Without the loss of time, measured pH data can be transferred directly to the Sleeplab System.

The Orion II SL is provided with an additional Lemo-3-Connector that can be connected to up to two analogue entries of the Sleeplab System.

The product range includes a 3 meter connection cable. The end of the cable has no plug and can be adapted to the connector of any Sleeplab System.

Warning
isolated
input!



The Orion II is without doubt one of the most comprehensive user friendly pH recorders available to diagnose in Gastro-Esophageal Reflux Diseases (GERD) and non cardiac chest pain.

■ Bluetooth®

Wireless Communication up to 10 meter.
Easy Calibration of catheter via the PC.
Automatic Data Transfer while patient is still in waiting room.

■ VIP, Virtual Instructor Program™

Guides user to successfully complete procedures with minimal training.

■ 4 pH channels (Orion II SL: 2 pH channels) Glass and antimony.

■ Lightweight

Compact housing with comfortable carrying bag.



Medical Measurement Systems

P.O. Box 580

7500 AN Enschede

The Netherlands

T + 31 - 53 - 480 37 00

F + 31 - 53 - 480 37 01

E info@mmsinternational.com

I www.mmsinternational.com