

Evac 40 Consumables and Accessories

The portable electrical aspirator Evac 40 is used for the aspiration and collection of secretions, body fluids, tissue parts and lavage fluids during surgical, endoscopic or minimally invasive procedures.



MT-0194-2005

Benefits

Evac 40

The Evac 40 is a high-performance surgical aspirator for high demands.

The vacuum pump with membrane technology is designed for continuous operation and distinguishes itself by its quiet running properties. With the foot switch the aspiration can be easily started or stopped. Especially easy-to-clean surfaces and an integrated, large precision vacuum display round out the system. An overflow safeguard in the lid of the 4 L secretion jar and a protection jar ensure that secretion cannot enter the device. The OR aspirator is provided with two rails, and can thus be used with different rail-mounted products.

Benefits at a glance:

- For all surgical fields of application
- Powerful, yet quiet
- Maintenance-free diaphragm pump
- Easy operation and care
- Mobile device on large castors

Technical Data

| | |
|--------------------------------------|---|
| Device suction power | 40 L/min |
| Max. vacuum | approx. – 90 kPa ¹⁾ |
| Vacuum display | -1 to 0 bar -750 to 0 mmHg |
| Secretion jar | 40 L/min |
| Rated voltage 230 V / 110 V | 230 V, 50/60 Hz; 110 V, 50/60 Hz |
| Rated current 230 V / 110 V | approx. 0,45 A at 230 V~/ approx. 0,82 A at 110 V |
| Rated output 230 V / 110 V | approx. 100 W / approx. 90 W |
| Operating time | Continuous operation within the stated voltage ranges |
| Dimensions (without secretion flask) | 999 x 435 x 436 mm |
| Weight (without secretion jar) | approx. 18 kg |
| Max. sound pressure | 60 dB (A) |
| Classification IX | Class II a (according to Directive 93/42/EEC Annex) |
| UMDNS code | 10-208 |
| ¹⁾ 1 bar = 1 kPa x 100 | |

Ordering Information

| Evac 40 | Order number |
|---|--------------|
| basic unit with foot switch, 230 V | 2M86752 |
| basic unit with foot switch, 110 V | MP00560 |
| Accessories necessary for operation: | |
| Secretion jar 4 L, incl. basket and cap | 2M85442 |
| Rail holder | 2M86757 |
| Suction filter (set of 10) | MX37021 |
| Optional accessories: | |
| Spare jar 4 L incl. basket | 2M85443 |
| Receiving jar with rail holder | 2M86758 |
| Connecting hose holder | 2M85446 |
| Aspiration hose holder | 2M85239 |

Notes

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
Moislinger Allee 53–55
23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany

As of August 2015

Dräger Medical GmbH changes
to Drägerwerk AG & Co. KGaA

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Medical GmbH
Moislinger Allee 53–55
23558 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Médical S.A.S.
Parc de Haute Technologie
d'Antony 2
25, rue Georges Besse
92182 Antony Cedex, France
Tel +33 1 46 11 56 00
Fax +33 1 40 96 97 20
d1mfr-contact@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Medical GmbH
Branch Office
P.O. Box 505108
Dubai, United Arab Emirates
Tel +971 4 4294 600
Fax +971 4 4294 699
contactuae@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama Comercial
S. de R.L.
Complejo Business Park,
V tower, 10th floor
Panama City
Tel +507 377 9100
Fax +507 377 9130
contactcsa@draeger.com

REGION ASIA / PACIFIC

Draeger Medical
South East Asia Pte Ltd.
25 International Business Park
#04-27/29 German Centre
Singapore 609916, Singapore
Tel +65 6572 4388
Fax +65 6572 4399
asia.pacific@draeger.com

REGION NORTH AMERICA

Draeger Medical, Inc.
3135 Quarry Road
Telford, PA 18969-1042, USA
Tel +1 215 721 5400
Toll-free +1 800 437 2437
Fax +1 215 723 5935
info.usa@draeger.com

Locate your Regional
Sales Representative at:
www.draeger.com/contact



The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).