FAMED
Żywiec

Member of REINSBERG® GROUP

We have been manufacturing
medical equipment of the highest quality since 1947.



BEDSIDE LOCKERS AND TABLES

product catalogue

The products are made of materials with antibacterial properties.





Famed Żywiec is a leading European manufacturer of high quality medical equipment. For over 75 years we have been supplying medical institutions with thousands of units of the highest quality operating tables, hospital and delivery beds or treatment chairs. Our products are used by doctors in Poland and in over 120 countries around the world.

We are distinguished by our qualified staff, wide portfolio of products, modern machinery and unique technological solutions used in production. We were one of the first companies in the world to introduce carbon fibre in the production of operating tables. We also offer the world's safest medical imaging table.

Our mission is to make the work of medical staff easier and to increase comfort for patients during their stay in hospital. Innovative projects and the pursuit of reliability in our products directly contribute to improving the quality of medical care in Poland and worldwide.

President of Famed Żywiec

Marek Suczyk

FAMED
Żywiec

Member of REINSBERG® GROUP

BEDSIDE LOCKERS AND TABLES

TABLE OF CONTENTS

1. BEDSIDE LOCKERS

Famed SP.02.0
Double-sided steel locker for storing patient's belongings with sliding, adjustable top

Famed SP.02.2
Simple locker for storing
patient's belongings with a shelf for bottles.

Famed **SP.02.0**

Very durable locker for storing patient's belongings with extendable table top with angle adjustment. Adjusting the height of the tabletop and its angle makes it possible to adapt the position of the tabletop to the patient's position on the bed.

CONSTRUCTION AND FINISH

This model of hospital locker is made of steel profiles and galvanised sheet steel, which form an ideal combination with the plastic tops. The drawer has been fitted with a two-chamber plastic insert and the cabinet box with an additional shelf for patient's belongings. Ergonomic wheels make it easy to position the locker in the hospital room, and the intuitive locking mechanism in two wheels ensures stability and prevents accidental change of position.

ADJUSTMENTS

Doors mounted on both sides of the locker and a double-sided drawer allow the locker to be positioned on both the left and right sides of the bed. An additional pull-out side table top with height and angle adjustment allows for easy and quick adaptation to the patient's position on the bed.

APPLICATIONS

It can be used both in hospital rooms and in care homes or the patient's home. The drawer and lockable locker allows storage of patient belongings, and the adjustable table can be used for meals or as a stand for a laptop or notebook. The lowest shelf has ventilation holes and can be used to store shoes.



Width 510 mm | **Depth** 390 mm | **Locker height** 880 mm | **Tabletop tilt angle** +/- 30° **Length of side top** 580 mm | **Width of side table top** 340 mm **High of side table top** 765-1010 mm | **Weight** 24 kg

ADJUSTMENT OPTIONS:



HIGH ADJUSTMENT

The locker has an additional side table top with height and angle adjustment. The table top is very stable and its positions can be changed quickly and easily to suit the height of the bed and the patient's requirements. Simply pull upwards on the black handle to adjust the positions in a smooth way. When it is not in use, the top can be easily hidden back - the table top hides into the side of the locker.



TILT ANGLE ADJUSTMENT

In addition to height adjustment, the side table top is also angle-adjustable. Using the red button, we can easily change its positioning according to the patient's needs. In this way, it can be used for meals or as a stand for a laptop, notepad or reading book.



SIDE CHANGE

Doors mounted on two sides of the locker allow access to the items inside, regardless of the locker being placed on the right or left side of the bed. A similar solution was used in the design of the drawer, which can be pulled out from two sides and a special locking mechanism integrated into the design prevents it from accidentally falling out.


Button and lever to adjust table position.


Lift the lever and move the tabletop to the required position.


Lift the tilt adjustment lock button on the tabletop and move it to the correct position.

Famed **SP.02.2**

A simple hospital locker for storing patient belongings. Additional space for bottles or other items in the side wall of the cabinet provides even more storage. A large box with ventilation holes and an additional shelf, as well as a capacious drawer, will keep the patient's personal items in order.

CONSTRUCTION AND FINISHING

This model of hospital locker is made of steel profiles and galvanized sheet steel, which make an ideal combination with the plastic tabletop. We have equipped the drawer with a two-chamber ABS plastic insert and the cabinet box with an additional shelf for the patient's personal items. Ergonomic wheels make it easy to position the locker in the hospital room, and a locking mechanism in two wheels ensures stability and prevents accidental changes of position.

ADJUSTMENT OPTIONS

The doors mounted on both sides of the locker and the double-sided drawer allow the locker to be positioned either on the left or right side of the bed.

APPLICATIONS

It can be used in hospital rooms as well as in nursing homes or the patient's home. The drawer and lockable cabinet allow storage of patient essentials, and the ABS shelf in the outer wall of the cabinet can be used to store bottles or items that the patient needs to have at hand.



Width 510 mm | **Depth** 390 mm | **Height** 880 mm | **Weight** 24 kg

TECHNICAL DATA COMPARISON

Model	hight	width	depth/length	weight	materials	other parameters
SP.02.0	880 mm (locker) 765-1010 mm (tabletop)	390 mm	510 mm	24 kg	steel	side tabletop: 580x340 mm, 30°
SP.02.2	880 mm	390 mm	510 mm	24 kg	steel	-

CUSTOMISATION OPTIONS

You have the opportunity to personalize your cabinets by choosing the colour of their fronts from our colour palette.*

Standard locker colours:

RAL 5002

RAL 5024

RAL 6027

RAL 1015



ADDITIONAL ACCESSORIES

- Bottle basket - WT-06**
- Towel holder with hooks - WT-07**

* Please note that the colours shown in the images may differ minimally from the actual colours, which is due to the difference in brightness of the display screens and differences in the printing of the product catalogues.
** Additional accessories are available for the SP-02.0 and SP -02.2 cabinets.



pSilver™

All ABS elements were made using pSilver™ technology. The prefabricate contains silver nanoparticles, which we use in the process of plastics production, provides effective bacteriostatic protection. The use of pSilver™ technology in the places that are most often touched by patients limits the proliferation of dangerous bacteria, fungi or viruses. The pSilver™ technology also involves a process of plastic element processing providing a smooth surface, whose texture prevents the accumulation of impurities.

FAMED
Żywiec

Member of REINSBERG GROUP

Famed Żywiec Sp. z o.o.
ul. Fabryczna 1
34-300 Żywiec, Poland

Office:
tel.: +48 33 866 62 00
fax: +48 33 475 58 90
sekretariat@famed.com.pl

Export:
tel.: +48 33 866 63 75
export@famed.com.pl

www.famed.com.pl



TÜVRheinland CERTIFIED
Management System
ISO 9001:2015
EN ISO 13485:2016
+AC:2018+A11:2021
www.tuv.com
ID 9000037697

EU Made in POLAND

UNEP/IE CLEANER PRODUCTION FRAMEWORK
CZYSTSZA PRODUKCJA

Edition 274/2026/10/5. Famed Żywiec reserves the right to modify the product and specification as part of technical progress. All illustrations and photos used in this material are for illustrative purposes only and may not reflect the finished product. The people shown in the photos are not medical professionals. They are models. The device presented in the catalog is intended for use in health care facilities by authorized persons after having read the manual (IFU).